

BRODA[®]

Hygiene Seating Solutions

Safety & Comfort for Hygiene Care



BRODA[®]

When it comes to hygiene care, the right shower commode chair can make the difference between a negative bathing experience and a positive one for both patient and caregiver. That's why Broda shower commodes have unique features to make hygiene comfortable and safe for everyone involved.

Table of Contents

Broda Benefits	1
Why Broda	3
Sentinel Shower Commode.....	4
Revive Shower Commode.....	5
Revive Bariatric Shower Commode.....	6
Configure Your Shower Commode.....	7
Accessories	7



Broda Benefits

At Broda, we promise to promote and preserve the wellbeing of the people who use our products. You can be sure your Broda chair or wheelchair will provide you with the following benefits.



Seating Comfort

We believe comfort is essential to a high quality of life and overall wellbeing.



Positioning

Our products offer postural alignment, pressure relief, and secure seating.



Mobility

Broda will help you reach the independence and social connection you deserve.



Patient Safety

Each of our products protects the user with built-in safety features



Caregiving Safety

We are dedicated to making healthcare facilities safe for every caregiver.



Infection Control

Keeping you healthy is our priority. Our products are easy to clean and disinfect.



Durability

Our long-lasting products are guaranteed by our unbeatable 10-year frame warranty.



Personalization

We'll build a wheelchair to match your style and provide the perfect fit.



Easy Order


Our team of seating and positioning experts are here to guide you through the process.

Lowering Stress in Hygiene Care

According to studies on caregiver burden, bathing assistance has been shown to be the most taxing aspect of care in nursing careers.¹

The most stressful aspects of bathing include...

- Up to five manual transfers
- Challenges dressing & undressing
- Awkward postures and movements
- Struggle to maintain upright patient posture



Hygiene Care Challenges

- Transferring a patient in order to bathe them may be difficult and uncomfortable
- Basic stock equipment may be uncomfortable, overused, worn out, and easily broken
- Patients & caregivers face fall risks in tight, slippery spaces
- Patients are vulnerable and need to maintain their dignity



1: Yamamoto, A., Azuma, T., Barrameda, H. C., Kida, N., Goto, A., & Ota, T. (2015). Caregiver's Eye Gaze and Field of View Presumption Method During Bathing Care in Elderly Facility. In V. G. Duffy (Ed.), Digital Human Modeling. Applications in Health, Safety, Ergonomics and Risk Management: Ergonomics and Health (Vol. 9185, pp. 524-532). Springer International Publishing. <https://doi.org/10.1007/978-3-319-21070-4>

The Broda Hygiene Difference

Safer Lifts and Transfers

Broda shower commodes enable safe transition options. Plus, the ease of maneuvering the patient reduces the need for adjustments and transfers.

Cleanliness Leads to Happiness

Comfort and dignity during the bathing process increases satisfaction with personal hygiene. This leads to better overall health and quality of life.

Cleanability

Easy-to-remove, minimally textured surfaces make Broda products and their components, such as bed pans and lateral supports, easy to clean.

Ergonomic Design

Broda shower commode chairs and accessories can adapt for a patient in multiple environments or stages of mobility, providing safety throughout the care continuum.

Durable Surfaces

Broda shower commodes feature durable, rust-resistant components that are easy to clean and disinfect between bathing sessions. Unlike basic shower commodes, Broda products provide safe, comfortable bathing solutions that last for years.

10-Year Frame Warranty

Quality products come with quality warranties. Broda guarantees our shower commode frames with an industry-leading 10-year warranty. While other shower commodes break frequently and need replacement, Broda shower commodes stand the test of time.



Ready for more information?

Scan here to view our shower/commode solutions.

How Broda Meets Hygiene Needs

- Dynamic tilt & recline provides safe, thorough, and efficient bathing
- Maximum durability with stainless steel frames and rust-resistant parts
- Removable arm supports create more space for safe transfers
- Unparalleled seating comfort reduces agitation and resistance to bathing

“ Broda’s Revive Shower Commode meets our patient’s needs. He indicates he feels safe and he is able to participate in his own care each morning. He is safe at all times while in the chair, and the caregivers are safe at all times while providing his care. All staff report being pleased with the commode and find it much easier to assist him. ”

—Pam
Director of Nursing





Sentinel

Shower Commode Wheelchair

RECOMMENDED FOR:



SPECIFICATIONS

BACK ANGLE	Angle adjustable to 5°, 10°, 15°
WEIGHT CAPACITY	Up to 300 lbs (136 kg)
SEAT WIDTHS	18", 20", 22" (46 cm, 51 cm, 56 cm)
SEAT DEPTHS	16", 17", 18", 19" (41 cm, 44 cm, 46 cm, 49 cm)
BACK HEIGHT	16", 18", 20" (41 cm, 46 cm, 51 cm)
SEAT SURFACE HEIGHTS	Adjustable 19.5" & 23.5" (50 cm, 60 cm)
WHEELS - MAG	24" Rear Mount (61 cm)
WHEELS - CASTER	5" (13 cm) Front total lock 5" (13 cm) Rear directional lock 5" (13 cm) Rear total lock
LIFETIME WARRANTY FRAME	16-Gauge Stainless Steel Frame 10 YEARS



All wheels can convert to casters




2-in-1 solid seat and commode



Configurable with either anti-rust, PVC mag wheels or cost-effective casters



Fold-up arm supports and swing-away leg supports



Scan to see accessories, order forms, and more!



Revive

Shower Commode Wheelchair

RECOMMENDED FOR:



SPECIFICATIONS

DEGREES OF TILT	Up to 12° (of posterior tilt) Up to 10° (of anterior tilt)
DEGREES OF RECLINE	Up to 67°
WEIGHT CAPACITY	Up to 400 lbs (181 kg)
SEAT WIDTHS	18", 22" (46 cm, 56 cm)
SEAT DEPTHS	18" (46 cm)
SEAT SURFACE HEIGHT	24.5" (63 cm)
BACK HEIGHT	26" (66 cm)
LIFETIME WARRANTY FRAME	16-Gauge Stainless Steel Frame 10 YEARS



Dynamic Tilting and Reclining



Removable Arm Supports for Easy Transfers



Anterior and Posterior Tilt



Rust-Resistant Frame and Components



Scan to see accessories, color options, order forms, and more!



Revive *Bariatric* Shower Commode Wheelchair

RECOMMENDED FOR:



SPECIFICATIONS

DEGREES OF TILT	Up to 14° (of posterior tilt) Up to 8° (of anterior tilt)
DEGREES OF RECLINE	Up to 65°
WEIGHT CAPACITY	Up to 800 lbs (363 kg)
SEAT WIDTHS	27", 33" (68 cm, 83 cm)
SEAT DEPTHS	21", 23" (54 cm, 58 cm)
SEAT SURFACE HEIGHTS	21.5", 23.5" (55 cm, 60 cm)
BACK HEIGHT	26" (66 cm)
LIFETIME WARRANTY FRAME	16-Gauge Stainless Steel Frame 10 YEARS



Concave back option accommodates and supports gluteal shelf




Dynamic tilting and reclining permits safe and thorough bathing



Anterior and posterior tilt allows for easy user positioning



Rust-resistant frame and components are durable and long-lasting




Scan to see accessories, color options, order forms, and more!

Configure Your Hygiene Wheelchair

Features, Functionality, and Customization Level

MORE INDEPENDENT



Weight Capacity: 300 lbs.
Tilt: 5°, 10°, 15°

MORE DEPENDENT



Weight Capacity: 400 lbs.
Posterior Tilt: Up to 12°
Recline: Up to 67° Recline

Bariatric



Weight Capacity: 800 lbs.
Posterior Tilt: Up to 14°
Anterior Tilt: Up to 8°
Recline: Up to 65°

Options & Accessories



Molded Lateral Supports



Commode pan with Basket



Pelvic Positioning Strap



Encore Rehab Wheelchair

Companion Wheelchair for Daily Use

The Encore Rehab Wheelchair offers a universal seating system, customizable for each user's seating and positioning needs. With 33 degrees of infinitely-adjustable tilt and a dynamic rocking feature, the Encore Rehab is perfect for a wide variety of care environments.



Scan to see accessories, color options, order forms, and more!



BRODA[®]

Broda USA

15455 Conway Rd, Suite 200
Chesterfield, MO 63017, USA

Broda Canada

560 Bingham Centre Dr
Kitchener, ON N2B 3X9, Canada

brodaseating.com

1-844-552-7632



Contract#
36F79721D0039

