


BRODA[®]

Parts & Accessories



CONTENTS

Anatomy of a Wheelchair	3
Lower Leg Supports	4
Arm Supports	7
Neck & Head Supports	8
Trays	9
Lateral Shoulder Support Cushions	10
Lateral Support Options	11
Positioning Straps	12
Thigh Belts	13
Wheels	14
Stabilizers and Anti-Tip	16
Accessories	17
Shower Commode Options	18
Cushions & Covers	20
Rehab Seating	30
Padding Packages	32
Kits	34
Cylinders & Cables	37

ANATOMY OF A WHEELCHAIR



Lower Leg Supports



Swing-Away Lower Leg Supports Black (Pair)

F-SW

500309050

Compatible with: Broda Pedal & Rehab Wheelchairs with Pin & Saddle Hangers



Swing-Away Lower Leg Supports with Angle-Adjustable Foot Support Plates Black (Pair)

F-SW-ANGLE

500309060

Compatible with: Broda Pedal and Rehab Wheelchairs with Pin & Saddle Hangers



Elevating Swing-Away Lower Leg Supports with Angle-Adjustable Foot Support Plates Black (Pair)

F-SW-ELV-ANGLE

500309065

Compatible with: Broda Pedal & Rehab Wheelchairs with Pin & Saddle Hangers



Elevating Swing-Away Lower Leg Supports Black (Pair)

F-SW-ELV

500309075

Compatible with: Broda Pedal & Rehab Wheelchairs with Pin & Saddle Hangers



Swing-Away Lower Leg Supports Silver (Pair)

F-SW-IKN

500309053

Compatible with: Encore Pedal Wheelchair, Encore Rehab Wheelchair, & Comfort Rehab Wheelchair with Drop-in Hangers



Elevating Swing-Away Lower Leg Support Encore Silver (Pair)

F-SW-ELV-IKN

500309078

F-SW-ELV-IKN-SML (16" seat width)

500309079

Compatible with: Encore Pedal Wheelchair, Encore Rehab Wheelchair, & Comfort Rehab Wheelchair with Drop-in Hangers



Swing-Away Lower Leg Support Encore Black with Angle-Adjustable Foot Support Plates Black (Pair)

F-SW-IKN-ADJ

500309090

F-SW-IKN-ADJ-SML (16" seat width)

500309091

Compatible with: Encore Pedal Wheelchair, Encore Rehab Wheelchair, & Comfort Rehab Wheelchair with Drop-in Hangers

Lower Leg Supports



**Swing-Away Lower Leg Support
75° Fixed Foot Support Plates
Black (Pair)**

*F-SW 75-587-V3
510009079*

Compatible with: Comfort Rehab
Wheelchair with Drop-In Hangers



**Swing-Away Lower Leg Support
60° Angle-Adjustable Foot
Support Plates Black (Pair)**

*F-SW 60ADJ-587V3
510009085*

Compatible with: Comfort Rehab
Wheelchair with Drop-In Hangers



**Swing-Away Lower Leg Support
75° Angle-Adjustable Foot
Support Plates Black (Pair)**

*F-SW-75ADJ-587V3
510009080*

Compatible with: Comfort Rehab
Wheelchair with Drop-In Hangers



**Elevating Swing-Away Lower
Leg Supports with Fixed Foot
Support Plates for 587V3 Black
(Pair)**

*F-SW-ELV-587V3
510009093*

Compatible with: Comfort Rehab
Wheelchair with Drop-In Hangers



**Elevating Swing-Away Lower Leg
Supports with Angle-Adjustable
Foot Support Plates 587V3 Black
(Pair)**

*F-SW-ELV-ANGLE-587V3
510009090*

Compatible with: Comfort Rehab
Wheelchair with Drop-In Hangers

Lower Leg Supports



Angle-Adjustable, Articulating Single Lower Leg Support

F-SINGLE

OptionID varies by width, strap color, and strap quantity

Compatible with: Broda Positioning Wheelchairs



Angle-Adjustable, Articulating Split Lower Leg Support

F-SPLIT

OptionID varies by width, strap color, and strap quantity

Compatible with: Broda Positioning Wheelchairs



Angle-Adjustable, Articulating Single Lower Leg Support CS385

F-SINGLE-CS385

OptionID varies by width, strap color, and strap quantity

Compatible with: Revive Shower Commode



Angle-Adjustable, Articulating Single Lower Leg Support 385B

F-SINGLE-385B

500309055 (grey straps)

500309056 (blue straps)

500309057 (green straps)

Compatible with: Revive Bariatric Shower Commode



Angle- and Length-Adjustable, Articulating ABS Lower Leg Support ABS

F-ABS-LEGFOOT

OptionID varies by width and frame color

Compatible with: Broda Positioning Wheelchairs

Factory install recommended.



Contracture Sling Frame 30VT

C-SLING-FRM-BC

OptionID varies by width and upholstery color

Compatible with: Centric Positioning Wheelchair



Contracture Sling Frame 785, 85V

C-SLING-FRM-ELV2/V3

OptionID varies by width and upholstery color

Compatible with: Elite Tilt Wheelchair, Elite Positioning Wheelchair




Lower Leg Support Leg and Foot Board 48, 48V4

LEG-FOOT-48/48R

OptionID varies by width

Compatible with: Pedal Wheelchairs, Latitude Rehab Wheelchair, Encore Rehab Wheelchair, & Comfort Rehab Wheelchair



Lower Leg Support Leg Board LEGPANEL

OptionID varies by width

Compatible with: Pedal Wheelchairs, Latitude Rehab Wheelchair, Encore Rehab Wheelchair, & Comfort Rehab Wheelchair

Arm Supports



Arm Support, Black (Each)

AR

570402000

Included in: Standard Chair Configurations (except Comfort Rehab Wheelchair)



Cushioned Arm Supports (Pair)

AR-PD

OptionID varies by color

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs & Tranquille Glider



Long Arm Support, Black (Pair)

AR-LONG

570402007

Compatible with: Synthesis Positioning, Rehab, & Transport Wheelchairs, Elite Positioning, Tilt, & Transport Wheelchairs



Arm Support Long Padded (Pair)

AR-LONG-PD

OptionID varies by color

Compatible with: Encore Rehab Wheelchair, Synthesis Positioning, Rehab & Transport Wheelchairs, Elite Positioning, Tilt & Transport Wheelchairs



Narrow Arm Support, Black (Each)

AR-N

570402001

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs & Tranquille Glider
Included in: Comfort Rehab Wheelchair



Arm Support Trough Triptex (Each)

AR-TR

100370075

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs & Tranquille Glider



Arm Support Fleece Cover (Pair)

AR-CVR-F

100100000

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs & Tranquille Glider

Neck & Head Supports




Neck Support

NR

OptionID varies by width and color

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs

Included in: Broda Standard Cushions



Neck Support Contoured Visco Cushion

NR-CV

OptionID varies by width and color

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs with seat width 16"-22"



Neck Support Soft Visco Cushion

NS

100270000

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs



Round Bolster Cushion - Each (pictured with Neck Support)

RB

OptionID varies by color

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs



Neck Support 100

NR-AL

OptionID varies by color

Compatible with: Tranquille Glider (high back option)



Molded 2" Offset Head Support CS385

HDRST-2"OFFSET

500800114

Compatible with: Revive Shower Commode & Revive Bariatric Shower Commode

Trays



Tray ABS with Strap 16" - 22"

TRAY-ABS

OptionID varies by tray color, width, and strap style

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs, Tranquille & large Aspire Gliders



Tray ABS with Strap 24" - 30"

TRAY-ABS-OVR

504300100 (24")

506300100 (26")

507300100 (28")

508300100 (30")

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs, Tranquille & large Aspire Gliders



Tray Clear with Strap 16" - 20"

TRAY-CLEAR

OptionID varies by width and strap style

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs, Tranquille & large Aspire Gliders



Tray ABS Cushioned Cover

TRAY-COVER

OptionID varies by width and color

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs, Tranquille & large Aspire Gliders




Half Tray with Lower Side Cover

TRAY-HALF

OptionID varies by tray color, upholstery color and left/right alignment

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs (except Centric Positioning Wheelchair), Tranquille & large Aspire Gliders

Replacement Tray Parts



Tray Slider Rails

TRAYSLIDERS

430306100

Compatible with: Oversize Trays 24"-30" widths



Tray 1" Strap

TRAY-STRAP

100800020 (center buckle)

100800022 (side buckle)

Compatible with: Broda Positioning, Pedal & Rehab Wheelchairs

Lateral Shoulder Support Cushions



Lateral Shoulder Support Cushion (Pair)

SB-SH

OptionID varies by color

Compatible with: Sashay Pedal Wheelchair, Latitude Pedal and Rehab Wheelchairs, Elite Wheelchairs, Midline & Centric Positioning Wheelchairs



Lateral Shoulder Supports V4 Cushion (Pair)

SB-SH-V4

OptionID varies by color

Compatible with: Synthesis and Encore Wheelchairs



Lateral Shoulder Supports 985 Cushion (Pair)

SB-SH-985

100301035 (Burgundy)

100302035 (Sapphire)

100307035 (Charcoal)

Compatible with: Vanguard Positioning Wheelchair



Long Lateral Shoulder Supports (Pair)

SB-LONG

OptionID varies by color

Compatible with: Sashay Pedal Wheelchair, Latitude Pedal and Rehab Wheelchairs, Elite Wheelchairs, Midline & Centric Positioning Wheelchairs



Long Lateral Shoulder Supports V4 (Pair)

SB-LONG-V4

OptionID varies by color

Compatible with: Synthesis and Encore Wheelchairs



Acute Care Lateral Shoulder Supports (Pair)

SB-AC-SH

OptionID varies by color

Compatible with: Sashay Pedal Wheelchair, Latitude Pedal and Rehab Wheelchairs, Elite Wheelchairs, Midline & Centric Positioning Wheelchairs



Long Acute Care Lateral Shoulder Supports (Pair)

SB-AC-LONG

OptionID varies by color

Compatible with: Sashay Pedal Wheelchair, Latitude Pedal and Rehab Wheelchairs, Elite Wheelchairs, Midline & Centric Positioning Wheelchairs



Round Lateral Shoulder Supports (Pair)

SB-RD-SH

OptionID varies by color

Compatible with: Sashay Pedal Wheelchair, Latitude Pedal and Rehab Wheelchairs, Elite Wheelchairs, Midline & Centric Positioning Wheelchairs with WING-ROD frame



Long Rounded Lateral Shoulder Supports (Pair)

SB-RD-LONG

OptionID varies by color

Compatible with: Sashay Pedal Wheelchair, Latitude Pedal and Rehab Wheelchairs, Elite Wheelchairs, Midline & Centric Positioning Wheelchairs with WING-ROD frame

Lateral Support Options



Removable Lateral Shoulder Supports Long Cover (Pair)
WC-LONG
100301060



Removable Lateral Shoulder Supports Short Cover (Pair)
WC-SH
OptionID dependent on upholstery color

Lateral Torso Supports



Molded Lateral CS385 Lateral Support (Pair)
LATERAL-MD
500800115
Compatible with: Revive Shower Commode, Revive Bariatric Shower Commode, Sentinel Shower Commode



Lateral Torso Supports (Pair)
LTS
500001201 (*grey*)
500001202 (*black*)
Compatible with: Positioning, Pedal, & Rehab Wheelchairs
Factory Installation Recommended



Comfort Company Lateral Torso Supports (Each)
CC-LTS
430400004
Compatible with: Positioning, Pedal & Rehab Wheelchairs with 16"-22" seat width

Positioning Straps



Pelvic Positioning Strap Black S-BELT

OptionID varies by size

Compatible with: Positioning,
Pedal & Rehab Wheelchairs &
Tranquille Glider



Pelvic Positioning Strap Black Auto Buckle

S-BELT-A

OptionID varies by size and buckle
Compatible with: Positioning,
Pedal & Rehab Wheelchairs &
Tranquille Glider



Pelvic Positioning Strap Cushion S-BELT-PD

OptionID varies by size and color

Compatible with: Positioning,
Pedal & Rehab Wheelchairs &
Tranquille Glider



Pelvic Positioning Strap Cushion Auto Buckle

S-BELT-PD-A

OptionID varies by size and color

Compatible with: Positioning,
Pedal & Rehab Wheelchairs &
Tranquille Glider



Pelvic Positioning Strap Black WC19

S-BELT-WC19

500001185

Included in: WC19 Transport
Package



2-Point Pelvic Belt

S-BELT-587

520001000

Compatible with: Comfort Rehab
Wheelchair

Strap Sizes

Small: 16", 18", 20" seat widths

Medium: 22", 24" seat widths

Large: 26", 28" seat widths

Thigh Belts



Thigh Belt Cushioned

T-BELT

OptionID varies by size and color

Compatible with: All Broda Wheelchairs (Except Hygiene & Gliders)

Included in: Huntington's Special Padding



Thigh Belt Cushioned Auto Buckle

T-BELT-A

OptionID varies by size and color

Compatible with: All Broda Wheelchairs (Except Hygiene & Gliders)



Thigh Belt Cushioned Neoprene

T-BELT-NEO

100550020 (small)

100550021 (medium)

100550022 (large)

Compatible with: All Broda Wheelchairs (Except Hygiene & Gliders)



Thigh Belt Cushioned Neoprene Auto Buckle

T-BELT-NEO-A

100500160 (small)

100500161 (medium)

100500162 (large)

Compatible with: All Broda Wheelchairs (Except Hygiene & Gliders)



Thigh Belt Cushioned Neoprene for APP Package

T-BELT-NEO-APP

100550023 (small)

100550024 (medium)

100550025 (large)

Included in: APP Package



Thigh Belt Cushioned Neoprene Auto Buckle for APP Package

T-BELT-NEO-A-APP

100500163 (small)

100500164 (medium)

100500165 (large)

Compatible with: APP Package

Strap Sizes

Small: 16", 18", 20" seat widths

Medium: 22", 24" seat widths

Large: 26", 28" seat widths

Wheels



Caster Wheel 5" Total Lock (Each)

CASTER-5"TL

500500001 (grey)

500500006 (black)

Available on: Broda Positioning & Pedal Wheelchairs



Caster Wheel 5" Directional Lock (Each)

CASTER-5"DL

500500003 (grey)

500500007 (black)

Available on: Centric Positioning Wheelchair, Midline Positioning Wheelchair



Wheel 5" Aluminum Non-Locking Front (Each)

WHEEL-5"ALUM

500500018

Available on: Broda Pedal, Rehab, & Transport Wheelchairs



Caster Wheel 5" Twin Locking Wheel (Each)

CASTER-5"TL-TW

500500015 (grey)

500500016 (black)

Available on: Elite & Synthesis Transport Wheelchairs



Caster Wheel 6" Multi Directional Lock (Each)

CASTER-6"DL

500500009 (grey)

500500014 (black)

Available on: Elite & Synthesis Models



Caster Wheel 5" Fixed Wheel (Each)

CASTER-5"FX

500500005 (grey)

500500008 (black)

Available on: Encore, Latitude, & Sashay Models



Caster Wheel 5" Total Lock Stainless Steel (Each)

CASTER-5"TL-SS

500500002

Available on: Sentinel, Revive & Revive Bariatric Shower Commodes



Caster Wheel 5" Stainless Steel Directional Lock (Each)

CASTER-5"DL-SS

500500004

Available on: Revive Shower Commode



Caster Wheel 5" Total Lock Stainless Steel 385B (Each)

CASTER-5"TLSS-B

430106001

Available on: Revive Bariatric Shower Commode

Wheels



Caster Wheel 8" Multi Directional Lock 985 Wheel (Each)

CASTER-8"DL-985VG
430105010

Compatible with: Vanguard Positioning Wheelchair



Caster Wheel 6" Total Lock 985 (Each)

CASTER-6"TL
430101010

Compatible with: Vanguard Positioning Wheelchair



Wheel for 587V3 (Each)

WHEEL-587V3
430102513 (5")
430102514 (6")
430102515 (8")

Compatible with: Comfort Rehab Wheelchair



Mag Wheel 5 Spoke with Sleeve and Hub Cap (Each)

MAGWHEEL5
500500019 (20")
500500020 (22")
500500021 (24")

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs



Mag Wheel 14" (Each)

MAGWHEEL5
500500022
Compatible with: Elite Wheelchairs, Encore Pedal Wheelchair, & Synthesis Wheelchairs



Wheel Transport 3" Glider (Each)

WHEEL3"
430101201

Compatible with: Aspire Glider, Tranquille Glider

Brake Components



Parking Brake with Clamp

WC-BRAKE
500500500
Included in: Mag Wheel Kits



Parking Brake Extension Handles Black

WC-BRAKE-EXT-HNDL
500500511
Optional for: Chairs with Mag Wheels

Stabilizers and Anti-Tip



Anti Tip Devices V4, 550SR, V4RW 85V, 785

ATIP-BAR-ELV

500700112 (gray)

500700114 (black)

Compatible with: Elite & Synthesis Models



Rear 30VT, 48, 48R, MID, 500RW, 500SR Anti Tipper

ATIP-BAR-REAR

500700110 (gray)

500700121 (black)

Compatible with: Encore, Latitude, Sashay, Centric, & Midline Models



Front 30VT, 48, 48R, MID, 500RW, 500SR Anti Tipper

ATIP-BAR-FTM

500700111 (gray)

500700122 (black)

Compatible with: Encore, Latitude, Sashay, Centric, & Midline Models



Anti Tip Devices 587V3

ATIP-BAR-587V3

530001116

Compatible with: Comfort Rehab Wheelchair



Rear Stabilizer Frame V4 & 85V

STABIL-R-ELV2/V3

500700101 (gray)

500700102 (black)

Compatible with: Elite & Synthesis Models



**Front/Rear Stabilizer Frame
Select Models and Elite Front Stabilizer Frame**

STABIL-FR

500700106 (gray)

500700107 (black)

Compatible with: select wheelchair configurations



Stabilizer Frame Lateral

STABIL-LAT

OptionID dependent on seat width and frame color

Compatible with: select wheelchair configurations

Accessories



Frame Oxygen Holder and Assembly

O2-HOLDER

OptionID dependent on chair model, color, and wheels




Frame Oxygen Holder Adapter Front and Rear

O2-MOUNT-F/R

500800107 (grey)

500800137 (black)

OptionID dependent on chair model, color, and wheels



Frame Oxygen, IV Holder Mount V4, V4RW

O2/IV-MOUNT-V4

500800138

Compatible with: Synthesis Positioning & Rehab Wheelchairs



IV Pole and Assembly

IV-POLE

OptionID dependent on chair model, color, and wheels




Extra IV Pole Holder

IV-HOLDER

500800131 (black)

500800101 (grey)



O²/IV Side Adapter 85V 785 550SR

O2/IV-MOUNT-ELV2/V3

500800106 (grey)

500800136 (black)

Compatible with: Elite Models



Oxygen Holder, Soft Carry Bag

O2-HOLDER-BAG

100075106



Accessory Bag

ACCESSORY-BAG

570000120

Accessory Bag Personalization

BAG-INITIALS



Cup Holder

CUP-HOLDER

430003134

Shower Commode Options

Seat Options



Seat Commode Cushioned B385
SEAT-PD-385B
430300003
Compatible with: Revive Bariatric Shower Commode




Seat Commode Cushioned CS385
SEAT-PD-CS385
431300005 (18" seat width)
433300005 (22" seat width)
Compatible with: Revive Shower Commode




Seat Large Hole Cushioned CS385
SEAT-PD-PERI-LG-CS385
433301002
Compatible with: Revive Shower Commode



Seat Small Hole Cushioned CS385
SEAT-PD-PERI-SM-CS385
431301001 (18" seat width)
433301001 (22" seat width)
Compatible with: Revive Shower Commode



Seat Transport Cushioned CS385
SEAT-PD-TS-CS385
431300006 (18" seat width)
433300006 (22" seat width)
Compatible with: Revive Shower Commode



Seat Transport Cushioned B385
SEAT-PD-TS-385B
430300004
Compatible with: Revive Bariatric Shower Commode



Sentinel Seat Commode
SEAT-PD-CS313
501800187 (18" seat width)
509800187 (20" seat width)
503800187 (22" seat width)
Compatible with: Sentinel Shower Commode

Shower Commode Options

Bedpans



Cushioned Seat Bedpan CS385/ B385

BEDPAN-PD-V2

500800120

Compatible with: Revive Shower
Commode, Revive Bariatric
Shower Commode



Cushioned Seat Bedpan for Roll Over Basket V3

BEDPAN-PD-V3

500800121

Compatible with: Revive Shower
Commode, Revive Bariatric
Shower Commode



Bedpan Commode Aqua

BEDPAN-AQUA

430003201

Compatible with: Revive Shower
Commode




Basket Frame CS385

BASKET-CS385-PD-V3

200001605

Compatible with: Revive Shower
Commode




Basket Frame B385

BASKET-385B

200001602

Compatible with: Revive Bariatric
Shower Commode



Basket Frame Aqua CS385

BASKET-CS385-AQUA

200001601

Compatible with: Revive Shower
Commode



Sentinel Basket Frame

COMMUNE-BSKT-CS313

501800188 (18" seat width)

509800188 (20" seat width)

503800188 (22" seat width)

Compatible with: Sentinel Shower
Commode



Commode Pan with Lid

BED-PAN-CS313

500801188

Compatible with: Sentinel Shower
Commode

Cushions & Covers

Vanguard Positioning Wheelchair



1. Back Cushions & Covers

Back Cushion 985

BP-985

OptionID dependent on seat width

Included in: Broda Standard Cushions

Back Vinyl Cover 985

BC-V-985

OptionID dependent on seat width and upholstery color

Included in: Broda Standard Cushions

Back Cushion and Vinyl Removable Cover 985

BP-VC-985

OptionID dependent on seat width and upholstery color

Included in: Broda Standard Cushions

Back Full Vinyl Cushion and Cover 985

FVC-BP-985

OptionID dependent on back shape, seat width, back length, and upholstery color

Included in: Acute Care Configuration

2. Seat Cushions & Covers

Seat Cushion 985

SP-985

OptionID dependent on seat width

Included in: Broda Standard Cushions

Seat Cushion and Vinyl Removable Cover 985

SP-VC-985

OptionID dependent on seat width and upholstery color

Compatible with: Vanguard Positioning Wheelchair

Seat Vinyl Cover 985

SC-V-985

OptionID dependent on seat width and upholstery color

Included in: Broda Standard Cushions

Seat Full Vinyl Cushion and Cover 985

FVC-SP-985

OptionID dependent on seat width, seat depth, and upholstery color

Cushions & Covers

3. Lower Side Covers

Lower Base Mount Side Cover 985 (Pair)

LSP-BM-985

100701013 (Burgundy)

100702013 (Sapphire)

100707013 (Charcoal)

Included in: Broda Standard Cushions

Lower Seat Mount Side Cover 985 (Pair)

LSP-SM-985

100701014 (Burgundy)

100702014 (Sapphire)

100707014 (Charcoal)

Included in: Broda Standard Cushions

4. Upper Side Covers

Upper Side Cover 985 (Pair)

USP-ELV2/V3/985

100701041 (Burgundy)

100702041 (Sapphire)

100707041 (Charcoal)

Included in: Broda Standard Cushions

5. Terrycloth Options

Back Terry Cover 985

BC-T-985

OptionID dependent on seat width and terry color

Back Cushion and Terry Cover 985

BP-TC-985

OptionID dependent on seat width and terry color

Seat Terry Cover 985

SC-T-985

OptionID dependent on seat width and terry color

Seat Cushion and Terry Cover 985

SP-TC-985

OptionID dependent on seat width and terry color

6. Calf Cushion

Calf Support Cushion 985

CP-985

OptionID dependent on seat width and upholstery color

Included in: Broda Standard Cushions



Cushions & Covers

- 1 Back Full Vinyl Non-Removable Cushion and Cover 100**
BP-VC-GLIDER
OptionID dependent on chair model and upholstery color
Compatible with: Tranquille Glider
- 2 Seat Full Vinyl Non-Removable Cushion and Cover 100**
SP-VC-GLIDER
OptionID dependent on chair model and upholstery color
Compatible with: Tranquille Glider
- 3 Side Cover 100 (Pair)**
SP-AL-GLIDER
OptionID dependent on upholstery color
Compatible with: Tranquille Glider
- 4 Seat Cushion Full Vinyl Non-Removable Cushion & Cover 100P**
FVC-BP-PEDGLIDER
OptionID dependent on upholstery color and seat width
Compatible with: Aspire Glider

Tranquille Glider



Aspire Glider



Cushions & Covers

Synthesis Positioning Wheelchair

1 Back Full Vinyl Cushion and Cover V4

FVC-BP-V4

OptionID dependent on seat width and upholstery color

Compatible with: Synthesis Positioning Wheelchair

Seat Full Vinyl Cushion and Cover V4, 785

FVC-SP-785

OptionID dependent on seat width, seat depth and upholstery color

Compatible with: Synthesis Positioning Wheelchair, Elite Positioning Wheelchair

2 Seat Full Vinyl Cushion and Notch Cover V4, 785

FVC-SP-N-785

OptionID dependent on seat width, seat depth and upholstery color

Compatible with: Synthesis Positioning Wheelchair, Elite Positioning Wheelchair

3 Lower Side Cover V4 (Pair)

LSP-V4

OptionID dependent

Compatible with: Synthesis Positioning Wheelchair
Included in:

4 Upper Side Cover V4 (Pair)

USP-V4

OptionID dependent

Compatible with: Synthesis Positioning Wheelchair
Included in:



Cushions & Covers

1 Seat Full Vinyl Cushion and Cover MID w/ Slot
FVC-SP-MID
OptionID dependent on seat width, seat depth and upholstery color
Compatible with: Midline Positioning Wheelchair

2 Back Full Vinyl Cushion and Cover MID, 785
FVC-BP-MID/785
OptionID dependent on seat width and upholstery color
Compatible with: Midline Positioning Wheelchair, Elite Positioning Wheelchair

3 Lower Side Cover MID (Pair)
LSP-MID
100701090 (Burgundy)
100702090 (Sapphire)
100707090 (Charcoal)
Compatible with: Midline Positioning Wheelchair

4 Upper Side Cover MID (Pair)
USP-MID/48V2
100701080 (Burgundy)
100702080 (Sapphire)
100707080 (Charcoal)
Compatible with: Midline Positioning Wheelchair, Sashay Pedal Wheelchair

5 Lower Side Cover 30VT (Pair)
LSP-BC
100701010 (Burgundy)
100702010 (Sapphire)
100707010 (Charcoal)
Compatible with: Centric Positioning Wheelchair

6 Upper Side Cover 30VT (Pair)
USP-BC
100701020 (Burgundy)
100702020 (Sapphire)
100707020 (Charcoal)
Compatible with: Centric Positioning Wheelchair



**Midline
Positioning
Wheelchair**



**Centric
Positioning
Wheelchair**

Cushions & Covers

- 1 Seat Full Vinyl Cushion and Cover 48R**
FVC-SP-48R
OptionID dependent on seat width, seat depth, and upholstery cover
Compatible with: Latitude Pedal Wheelchair
- 2 Back Full Vinyl Cushion and Cover 48R**
FVC-BP-48R
OptionID dependent on seat width and upholstery color
Compatible with: Latitude Pedal Wheelchair
- 3 Lower Side Cover 48 V2 (Pair)**
LSP-48V2/48R
OptionID dependent on upholstery color
Compatible with: Sashay and Latitude Pedal Wheelchairs
- 4 Base Apron Cushion 48R**
APRON-48R
OptionID dependent on seat width, seat depth, and upholstery color
Compatible with: Latitude Pedal Wheelchair
- 5 Lower Side Cover 587V3 (Pair)**
LSP-587V3
OptionID dependent on upholstery color
Compatible with: Comfort Rehab Wheelchair

Latitude Pedal Wheelchair



Comfort Rehab Wheelchair



Cushions & Covers

1 Back Vinyl Cover Encore Short Back Cover

BC-V-ENCORE-SHORT

OptionID dependent on seat width and upholstery color

Compatible with: Encore Pedal Wheelchair

Back Terry Cover Encore Short Back

BC-T-ENCORE-SHORT

OptionID dependent on seat width and terry color

Compatible with: Encore Pedal Wheelchair

Back Full Vinyl Cushion and Cover Encore Short Back

FVC-BP-ENCORE-SHORT

OptionID dependent on seat width and upholstery color

Compatible with: Encore Pedal Wheelchair

Back Cushion Encore Short Back (Each)

BP-ENCORE-SHORT

OptionID dependent on seat width and upholstery color

Compatible with: Encore Pedal Wheelchair

2 Seat Full Vinyl Cushion and Cover Encore 500RW

FVC-SP-48RV4

OptionID dependent on seat width, seat depth, and upholstery color

Compatible with: Encore Pedal Wheelchair, Encore Rehab Wheelchair

Seat Full Vinyl Cushion and Notch Cover Encore 500RW

FVC-SP-N-48RV4

OptionID dependent on seat width, seat depth, and upholstery color

Compatible with: Encore Pedal Wheelchair, Encore Rehab Wheelchair

3 Back Full Vinyl Cushion and Cover Encore

FVC-BP-ENCORE

OptionID dependent on seat width and upholstery color

Compatible with: Encore Pedal Wheelchair

4 Seat Full Vinyl Cushion and Notch Cover Encore 48V4

FVC-SP-N-48V4

OptionID dependent on seat width and upholstery color

Compatible with: Encore Pedal Wheelchair

Seat Full Vinyl Cushion and Cover Encore 48V4

FVC-SP-48V4

OptionID dependent on seat width, seat depth, and upholstery color

Compatible with: Encore Pedal Wheelchair

5 Lower Side Cover Encore (Pair)

LSP-ENCORE

OptionID dependent on upholstery color

Compatible with: Encore Pedal Wheelchair

Encore Pedal Wheelchair



Cushions & Covers

1 Lower Seat Mount Side Cover 85V, 785 V2 (Pair)
LSP-SM-ELV2
100701033 (Burgundy)
100702033 (Sapphire)
100707033 (Charcoal)
Compatible with: Elite Tilt Wheelchair, Elite Positioning Wheelchair

Lower Base Mount Side Cover 85V, 785 (Pair)
LSP-BM-ELV2/V3
100701032 (Burgundy)
100702032 (Sapphire)
100707032 (Charcoal)
100720032 (Meteor)
Compatible with: Elite Tilt Wheelchair, Elite Positioning Wheelchair


Lower Seat Mount Side Cover 85V, 785 V3 (Pair)
LSP-SM-ELV3
100701034 (Burgundy)
100702034 (Sapphire)
100707034 (Charcoal)
100720034 (Meteor)
Compatible with: Elite Tilt Wheelchair, Elite Positioning Wheelchair

Elite Positioning Wheelchair



2 Upper Side Cover 85V, 785, 985 (Pair)
USP-ELV2/V3/985
100701041 (Burgundy)
100702041 (Sapphire)
100707041 (Charcoal)
100720041 (Meteor)
Compatible with: Elite Tilt Wheelchair, Elite Positioning Wheelchair, Vanguard Positioning Wheelchair


Cushions & Covers



Seat Vinyl Cover

SC-V

OptionID dependent on seat width, seat depth and upholstery color



Seat Cushion and Vinyl Removable Cover

SP-VC

OptionID dependent on seat width, seat depth and upholstery color



Seat Terry Cover

SC-T


OptionID dependent on seat width, seat depth, and terry color



Seat Cushion and Terry Cover

SP-TC

OptionID dependent on seat width, seat depth, and terry color



Back Vinyl Cover

BC-V

OptionID dependent on seat width and upholstery color



Back Cushion and Vinyl Removable Cover

BP-VC

OptionID dependent on seat width and upholstery color



Back Terry Cover

BC-T

OptionID dependent on seat width and terry color



Back Cushion and Terry Cover

BP-TC

OptionID dependent on seat width and terry color

Cushions & Covers



Side Reduction Cushion (Pair)

SR

OptionID dependent on upholstery color

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs



Side Arm Support Cushion (Pair)

SAS

OptionID dependent on upholstery color

Compatible with: Broda Positioning, Pedal, & Rehab Wheelchairs



Long Lower Leg Support Footbox with Pocket

FB-LONG-P

OptionID dependent on seat width and upholstery color

Compatible with: Broda Positioning Wheelchairs



Long Lower Leg Support Footbox

FB-LONG

OptionID dependent on seat width and upholstery color

Compatible with: Broda Positioning Wheelchairs
Included in: Huntington's Special Padding

Rehab Seating



Elite Back Support

ELITE-BACK-SUPPORT-22-26

ELITE-BACK-SUPPORT-16-20

MATRX-BACK-ELITE

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs



E2 Back Support

E2-BACK-SUPPORT-16-20

MATRX-BACK-E2

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs



PB Back Support

PB-BACK-SUPPORT-16-20

MATRX-BACK-PB

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs

Elite TR Back Support

ELITE-TR-BACK-SUPPORT-22-24

ELITE-TR-BACK-SUPPORT-16-20

MATRX-BACK-ELITE-TR

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs

E2 Deep Back Support

E2-BACK-SUPPORT-DEEP

MATRX-BACK-DEEP-E2

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs

PB Deep Back Support

PB-BACK-SUPPORT-DEEP-16-20

MATRX-BACK-DEEP-PB

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs

Elite Deep Back Support

ELITE-BACK-SUPPORT-DEEP-22-24

ELITE-BACK-SUPPORT-DEEP-16-20

MATRX-BACK-DEEP-ELITE

OptionID dependent on seat width and back height

Compatible with: Broda Rehab Wheelchairs

E2 Back Support Infection Control Slipcover

E2-SLIP-COVER-OFF

430309930

Compatible with: E2 Back Supports

E2 Deep Back Support Infection Control Slipcover

E2-DEEP-SLIP-COVER-OFF

430309932

Compatible with: E2 Deep Back Supports

Rehab Seating



Libra Seat Cushion

MATRIX-SEAT-LIBRA

OptionID dependent on seat width, seat depth and cushion cover

Compatible with: Broda Rehab Wheelchairs



Vi Seat Cushion

MATRIX-SEAT-VI

OptionID dependent on seat width, seat depth and cushion cover

Compatible with: Broda Rehab Wheelchairs



PSP Seat Cushion

MATRIX-SEAT-PSP

OptionID dependent on seat width, seat depth and cushion cover

Compatible with: Broda Rehab Wheelchairs

Libra Cushion Cover

MATRIX-COVER-LIBRA

OptionID dependent on seat width and depth

Compatible with: Libra Cushions

Vi Cushion Cover

MATRIX-COVER-VI

OptionID dependent on seat width and depth

Compatible with: Vi Cushions

PSP Cushion Cover

MATRIX-COVER-PSP

OptionID dependent on seat width and depth

Compatible with: PSP Cushions



PS Seat Cushion

MATRIX-SEAT-PS

OptionID dependent on seat width, seat depth and cushion cover

Compatible with: Broda Rehab Wheelchairs



Elan Head Support Hardware

MATRIX-ELAN-HRDW

Compatible with: Broda Rehab Wheelchairs



Loxx Head Support Hardware

MATRIX-LOXX-HRDW

Compatible with: Broda Rehab Wheelchairs

PS Cushion Cover

MATRIX-COVER-PS

OptionID dependent on seat width and depth

Compatible with: PS Cushions

Additional Matrix Options include:

Padding Packages

*Shown on Synthesis
Positioning Wheelchair*



Additional Positioning Padding

APP Package for Positioning Chairs

APP-POSITIONING-OFF
APP-POSITIONING-ON (on-chair)
OptionID dependent on seat width, seat depth, and back height
Compatible with: Positioning Wheelchairs

APP Package for Pedal Chairs

APP-PEDAL-OFF
APP-PEDAL-ON (on-chair)
OptionID dependent on seat width, seat depth, and back height
Compatible with: Pedal Wheelchairs

Additional Positioning Padding Includes:

APP Seat Cushion
APP Back Support Cushion Lower
APP Back Support Cushion Upper
APP Lateral Shoulder Supports
APP Head Support
APP Calf and Foot Support

Options available for individual purchase.

*Shown on Encore
Pedal Wheelchair*



Slim Positioning Padding

SPP Package

SPP-OFF
SPP-ON (on-chair)
OptionID dependent on seat width, seat depth, and back height
Compatible with: All Broda Wheelchairs except Shower Commodes and Gliders

Slim Positioning Padding Includes:

SPP Back Support Cushion Lower
SPP Back Support Cushion Upper (optional, not pictured)
SPP Seat Cushion

Options available for individual purchase.

Padding Packages

*Shown on Synthesis
Positioning Wheelchair*



Huntington's Special Padding

Huntington's Special Padding Replacement

HSP-PKG-UPGRD

OptionID dependent on chair model, seat width, and upholstery color

Compatible with: All Broda Wheelchairs except Shower Commodes, Gliders, Comfort Rehab Wheelchair, and Vanguard Positioning Wheelchair. Footbox not compatible on Encore models.

Huntington's Special Padding Includes:

Long Shoulder Supports
Cushioned Thigh Belt
Long Footbox

Options available for individual purchase.

*Shown on Synthesis
Positioning Wheelchair*



Acute Care Configuration

Acute Care Configuration (ACC) Package

ACC

OptionID dependent on chair model, seat width, and upholstery color

Compatible with: Synthesis Positioning, Rehab, & Transport Wheelchairs, Elite Positioning, Rehab & Tilt Wheelchairs, Midline Positioning Wheelchair, Centric Positioning Wheelchair, & Vanguard Positioning Wheelchair

Acute Care Configuration Includes:

Full Vinyl Wrap Seat and Back Cushions
ABS Lower Leg Support with Calf and Foot Cushions
ACC Lateral Shoulder Supports (except on Synthesis Models)

Options available for individual purchase.

Kits



Mag Wheel Conversion Kit for V4
MAG-CONV-V4
590005045 (20" wheel)
590005046 (22" wheel)
590005047 (14" wheel)
Compatible with: Synthesis
Positioning Wheelchair

Mag Wheel Conversion Kit for 48V4 & 48V4R
MAG-CONV-ENCORE
590005050 (20" wheel)
590005051 (22" wheel)
Compatible with: Encore Pedal
Wheelchair

Mag Wheel Conversion Kit for 785, 85V
MAG-CONV-ELV2/V3
*OptionID dependent on wheel
size and frame color*
Compatible with: Elite Models

Mag Wheel Conversion Kit for 48, 48R, 30VT, or MID
MAG-CONV
*OptionID dependent on wheel
size, placement, and frame color*
Compatible with: Latitude
Pedal Wheelchair, Sashay Pedal
Wheelchair, Centric Positioning
Wheelchair, Midline Positioning
Wheelchair

*Shown on Synthesis
Positioning Wheelchair*



Mag Wheel Conversion Kit for V4 - Middle Mount
MAG-CONV-V4-MIDDLEMNT
590005048
Compatible with: Synthesis
Wheelchairs

Mag Wheel Conversion Kits Include:
Mag Wheels
Brake
Brake Rod (where applicable)
Brackets
Anti-Tippers
Hardware

Middle Mount Mag Wheel Conversion for 785, 85V
MAG-CONV-ELV2/V3-MIDDLEMNT
590005034 (grey)
590005038 (black)
Compatible with: Elite Models

Options available for individual purchase.

*Shown on Synthesis
Positioning Wheelchair*

Kits



Back Frame Height Extension Kit

1 Back Support Height Kit Extension for All Models Kit

HGHT-BACK-KIT-3.0 (3" extension)

HGHT-BACK-KIT-1.5 (1.5" extension)

OptionID dependent on seat width, back height, and color combination

Compatible with: All Models

Seat Depth Frame Extension Kit

2 1.5" Seat Depth Extension Kit available on:

Elite Wheelchairs

Synthesis Wheelchairs

Midline Positioning Wheelchair

Centric Positioning Wheelchair

Latitude Pedal Wheelchair

Encore Pedal Wheelchair

OptionID and Alpha Code dependent on chair model, seat depth, seat width, and frame and upholstery color combination

3" Seat Depth Extension Kit available on:

Elite Wheelchairs

Midline Positioning Wheelchair

Centric Positioning Wheelchair

Latitude Pedal Wheelchair

Encore Pedal Wheelchair

OptionID and Alpha Code dependent on chair model, seat depth, seat width, and frame and upholstery color combination

Seat Surface Frame Height Kit

3 1.5" Seat Surface Height Kit for 85V, 785

HGHT-SEAT-KIT-EL-1.5

500300300 (grey)

500300310 (black)

Compatible with: Elite Models

3" Seat Surface Height Kit for 785, 85V

HGHT-SEAT-KIT-EL-3.0

500300305 (grey)

500300315 (black)

Compatible with: Elite Models

Seat Surface Height Kit for V4

HGHT-SEAT-KIT-V4

500300360 (19.5" seat height)


500300361 (20.5" seat height)

500300362 (22" seat height)

500300363 (23.5" seat height)

Compatible with: Synthesis Wheelchairs

Kits



Conversion Kit Lower Leg Support to Contracture 30VT

LGRT-CONV-C-30VT

OptionID dependent on upholstery color and seat width

Compatible with: Centric Positioning Wheelchair

Conversion Kit Contracture to Lower Leg Support 30VT

C-CONV-LGRT-30VT

OptionID dependent on seat width and upholstery color

Compatible with: Centric Positioning Wheelchair

Conversion Kit Lower Leg Support to Contracture 85V

LGRT-CONV-C-85V3

OptionID dependent on upholstery color and seat width

Compatible with: Elite Tilt Wheelchair

Cylinders & Cables



Cable and Handle Assembly
CABLE-HAND
OptionID dependent on length

Double Cable and Handle Assembly
CABLE-HAND-DBL
OptionID dependent on length

Cable and Handle Assembly 48R Pedal Rocker
CABLE-HAND-48R
 500600102
Compatible with: Latitude Pedal & Rehab Wheelchairs

Cable and Long Handle Assembly 48"
CABLE-HAND-LONG
 500600150

Cylinder Extension Bar
CYL-EXT-BAR
 500811004 (gray)
 500813004 (black)

Cable and Handle Assembly 587
CABLE-HAND-587
 500600152
Compatible with: Comfort Rehab Wheelchair

Double Cable and Handle with 32" Reversed Cable Length for PAB Assembly 985
CABLE-HAND-DBL-32RV-985
 500600107
Compatible with: Vanguard Positioning Wheelchair

Cable, Handle & Actuator Assembly
CABLE-ASSEMBLY-KIT
OptionID dependent on length and handle

Acuator and Lever
CABLE-ACT
 510006155

Cylinder and Actuator Kit

<i>CYL-LOCK-KIT-270N, PAB4</i>	500600116
<i>CYL-LOCK-KIT-300N, PAB5</i>	500600115
<i>CYL-LOCK-KIT-200N, PAB16</i>	500600146
<i>CYL-LOCK-KIT-600N, PAB17</i>	500600148
<i>CYL-LOCK-KIT-800N, PAB10</i>	500600144
<i>CYL-LOCK-KIT-800N-587, PAB6</i>	500600132
<i>CYL-LOCK-KIT-300N, PAB7, stainless steel</i>	500600128
<i>CYL-LOCK-KIT-600N, PAB8, stainless steel</i>	500600123
<i>CYL-LOCK-KIT-800N, PAB9, stainless steel</i>	500600125
<i>CYL-LOCK-KIT-270N, PAB11, stainless steel</i>	500600129
<i>CYL-LOCK-KIT-300N-385B, PAB12, stainless steel</i>	500600131
<i>CYL-LOCK-KIT-600N-48R, PAB13, stainless steel</i>	500600127
<i>CYL-HYD-985V2</i>	

