

PO: _____ **Date:** _____ **Customer #:** _____

Customer Name: _____ Ship To Name: _____
 Address: _____ Address: _____
 Phone: _____ Phone: _____
 Contact: _____ Contact: _____
 Email: _____ Tag: _____



Standard Features	
• Up to 14° Infinitely Adjustable Seat Tilt (Posterior)	• Up to 800lb Maximum Occupant Weight
• Up to 8° Infinitely Adjustable Seat Tilt (Anterior)	• Contoured Backrest
• Up to 65° Infinitely Adjustable Back Support Recline	• Comfort Tension Seating®
• Stainless Steel, Powder Coated Frame	• Width Adjusting Removable Arm Supports
• 4 x 6" Total Lock Stainless Steel Front & Rear Caster Wheels	
• 2 x 8" Locking Middle Wheels with One-Touch Emergency Brake	

Revive Bariatric Shower Commode - Model: B385 MSRP \$7,537.69
Preset Factory Installed Cylinder Lift Forces(Select one)

Cylinder Configuration	Suggested Occupant Weight (lbs)	Cylinder Lift Force	
2 x 300N	<input type="radio"/> 100-350 lbs.	135 lbs.	N/C
2 x 600N	<input type="radio"/> 350-500 lbs.	270 lbs.	N/C
2 x 800N	<input type="radio"/> 500-650 lbs.	360 lbs.	N/C
2 x 300N + 2 x 600N	<input type="radio"/> 650-800 lbs.	405 lbs.	N/C

*Guideline Only. Contact Broda for additional weight accommodation options

Seat and Back Support Options			
Seat Width:	<input checked="" type="radio"/> 27" N/C	Seat Depth:	<input type="radio"/> 21" N/C
	(33" with width adjustable armrests)		<input type="radio"/> 23" N/C
Seat Surface	<input type="radio"/> 21" N/C	Back Height:	<input checked="" type="radio"/> 26" Standard
Height:	<input type="radio"/> 23" N/C		
		No Anterior Tilt*	<input type="radio"/> N/C

*Max Upright Position - 5 degree posterior tilt

Wheel Options	
<input checked="" type="radio"/> 6" Total Lock Stainless Steel Front & Rear Caster Wheels & 8" Locking Middle Wheels with One-Touch Emergency Brake	Standard

Seat Options (Select one)	Bedpan & Basket Options
<input type="radio"/> Seat Commode Cushioned Standard	<input checked="" type="radio"/> Bedpan with Basket Standard
<input type="radio"/> Seat Transport Cushioned N/C	

Arm Support Options:	
Arm Support Options:	<input checked="" type="radio"/> Arm Supports Standard
Arm Support Style:	<input checked="" type="radio"/> Seat Width Extending (extends width to 33") Standard

Head Support	
Head Support:	<input checked="" type="radio"/> Molded 2" Offset Standard

Strap Color			
Strap Color:	<input type="radio"/> Grey	<input type="radio"/> Blue	<input type="radio"/> Green N/C

Common Accessories	
Lateral Supports:	<input type="radio"/> Molded Laterals (pair) \$92.00
Pelvic Belt:	<input type="radio"/> Pelvic Belt (side buckle)(large size) \$61.34

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
 - Residential Delivery Service:** Delivery to a residential location
 - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
 - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
 - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**