

PO: _____

Date: _____

Customer #: _____

Customer Name: _____

Ship To Name: _____

Address: _____

Address: _____

Phone: _____

Phone: _____

Contact: _____

Contact: _____

Email: _____

Tag: _____



Standard Features	
• 2 x 24" Rear Mag Wheels	• Swing Away Footrests w/ Heel Loops (14"-18" knee-to-heel range)
• 2 x 5" Total Lock Front Stainless Caster Wheels	• 10° Back Angle (adjustable to 5° and 15°)
• Adjustable Seat Surface Height	• Swing Up Arm Supports
• Adjustable Seat Depth	• Stainless Steel Adjustable Anti-tippers
• Adjustable Rear Mag Position (14"-21")	• 300lb Maximum Occupant Weight
• Adjustable Seat Angle (0°-6°)	• Stainless Steel, Powder Coated Frame
• Tension Adjustable Back Straps	

Sentinel Shower Commode - Model: CS313 MSRP \$2200.00

Seat and Back Support Options

Seat Width:	<input type="radio"/> 18" Standard	Seat Depth:	<input type="radio"/> 17" Standard
	<input type="radio"/> 20" \$400		<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 19" \$100
	<input type="radio"/> 22" \$500		
Seat Height:	<input type="radio"/> 19.5" Standard	Back Height:	<input type="radio"/> 18" Standard
	<input type="radio"/> 23.5" N/C		<input type="radio"/> 16" <input type="radio"/> 20" N/C

Parking Brake Extension Handles \$78.51

Front Options	Wheel Options	Rear Options
---------------	---------------	--------------

<input checked="" type="radio"/> 5" Front Total Lock Stainless Caster Wheels Standard	<input type="radio"/> 24" Mag Wheels- Pair Standard	
	2 rear caster options below: (\$400)	
	<input type="radio"/> 5" Total Lock Stainless Caster Wheels (2) N/C	
	<input type="radio"/> 5" Directional Lock Stainless Steel Caster Wheel (1) & 5" Total Lock Stainless Caster Wheel (1)	

Seat Options	Bedpan Options
--------------	----------------

<input checked="" type="radio"/> Two in One Commode/Shower Seat Standard	<input type="radio"/> Bedpan with Basket \$62
--	---

Common Accessories

Lateral Supports:	<input type="radio"/> Molded Laterals - Pair (pre-cut slots in back included)	\$92.00
Pelvic Belt:	<input type="radio"/> Pelvic Belt (18", 20" seat widths = small, 22" = medium)	\$61.34
	<input type="radio"/> Center Buckle <input type="radio"/> Side Buckle	
<input type="radio"/> Back Rest Pad with Slots 20":		\$50.00

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
 - Residential Delivery Service:** Delivery to a residential location
 - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
 - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
 - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours, Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**