

PO: \_\_\_\_\_

Date: \_\_\_\_\_

Customer #: \_\_\_\_\_

Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Email: \_\_\_\_\_

Ship To Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Tag: \_\_\_\_\_



Standard Features	
• Comfort Tension Seating®	• 250lb Maximum Occupant Weight
• Tubular Steel, Powder Coated Frame	• 2 x 20" Rear Mag Wheels
• Adjustable Seat Surface Height	• 2 x 5" Front Total Lock Caster Wheels
• Height Adjustable Arm Supports	• Up to 25° Infinitely Adjustable Seat Tilt
• Angle Adjustable, Articulating Lower Leg Support	• Up to 45° Infinitely Adjustable Back Recline

**Centric Positioning Wheelchair - Model: 30VT MSRP \$4,017.24**

Seat and Back Support Options			
<b>Seat Width:</b>	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28"	N/C \$484.53 \$1,097.84	
<b>Seat Surface Height:</b>	<input type="radio"/> 19.5" <input type="radio"/> 20.5" <input type="radio"/> 21.5"	N/C Standard N/C	
<b>Seat Depth:</b>	<input type="radio"/> 15.5" <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$116.54 Standard N/C \$116.54	
<b>Back Support Height:</b>	<input type="radio"/> 32" <input type="radio"/> 33.5" <input type="radio"/> 35"	Standard \$116.54 \$116.54	

Front Options	Wheel Options	Rear Options
<input type="radio"/> 2x 5" Total Lock Caster Wheels - Pair	Standard	<input type="radio"/> 20" Mag Wheels - Pair
<input type="radio"/> 20" Mag Wheels - Pair	N/C	<input type="radio"/> 1x 5" Directional Locking 1x Total Lock Caster Wheel
<input type="radio"/> 22" Mag Wheels - Pair	N/C	<input type="radio"/> 22" Mag Wheels - Pair
<input type="radio"/> 24" Mag Wheels - Pair	N/C	<input type="radio"/> 24" Mag Wheels - Pair
		<input type="radio"/> Parking Brake Extension Handles - Pair
		Standard (\$644.28) N/C N/C \$78.50

Lower Leg Support Options (Select One)			
<input type="radio"/> Angle Adjustable, Articulating Split Lower Leg Support	(not available in 16" seat width)		Standard
Choose Number of Straps (support length range)	<input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")		
<input type="radio"/> Angle Adjustable, Articulating Single Lower Leg Support			N/C
Choose Number of Straps (support length range)	<input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")		
<input type="radio"/> Contracture Sling (in place of lower leg support)			\$306.67

Arm Support Options			
<b>Arm Support:</b>	<input type="radio"/> Soft Armrest	Standard	<input checked="" type="radio"/> Height Adjustable Arm Supports
<b>Cushion:</b>	<input type="radio"/> Arm Trough (each)	\$159.46	<input checked="" type="radio"/> Base Mounted Arm Supports
<b>Options:</b>	<input type="radio"/> Cushioned Flip Arm Supports	\$142.05	<b>Mount:</b>
	<input type="radio"/> Cushioned Arm Supports* (pair)	\$92.00	
	*Not compatible with a tray		

Seating and Frame Color Options							
<b>Strap Color:</b>	<input type="radio"/> Grey	<input type="radio"/> Blue	N/C	<b>Cushion</b>	<input type="radio"/> Burgundy	<input type="radio"/> Charcoal	N/C
	<input type="radio"/> Green		N/C	<b>Color:</b>	<input type="radio"/> Sapphire Blue		N/C
<b>Frame Color:</b>	<input checked="" type="radio"/> Grey		Standard				

Cushions - Choose Option 1, 2, OR 3						
<input type="radio"/> <b>Option 1 - Broda Cushions</b>						N/C
Standard Broda Cushions:						
<b>Includes:</b> Neck Support, Lateral Shoulder Supports, Foot Support Cushions, Upper & Lower Side Covers, Seat & Back Cushions with Covers						
<b>Seat &amp; Back Cushion</b>	<input type="radio"/> Grey Terry Cover					N/C
<b>Cover Options:</b>	<input type="radio"/> Vinyl					N/C
	<input type="radio"/> Startex Cover					N/C
<b>Lateral Shoulder Support</b>	<input type="radio"/> Standard Short	<input type="radio"/> Rounded Short	<input type="radio"/> ACC Short			N/C
<b>Options:</b>	<input type="radio"/> Standard Long	<input type="radio"/> Rounded Long	<input type="radio"/> ACC Long			\$159.46
<b>Huntington's Special</b>	<input type="radio"/> Includes Long Lateral Shoulder Support Covers, Thigh Belt & Long Lower Leg Support Footbox					\$631.72
<b>Padding (HSP) package:</b>						

**Optional Substitutions Available Between APP & SPP Components**

<input type="radio"/> <b>Option 2 - Broda Additional Positioning Padding (APP)</b> (not compatible with 35" back support height)								
Broda APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Covers								
<b>Seat Width:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	<b>Seat Depth:</b>	<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$640.31
	<input type="radio"/> 22"	<input type="radio"/> 24"			<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$770.33
	<input type="radio"/> 26"	<input type="radio"/> 28"			<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$900.35

<input type="radio"/> <b>Option 3 - Broda Slim Posture Padding (SPP)</b>									
Broda SPP Padding Options: Seat, Back Support, and Back Support Extension Cushions									
<b>Seat Width:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	<b>Seat Depth:</b>	<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$211.98	
	<input type="radio"/> 22"	<input type="radio"/> 24"			<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$233.42	
	<input type="radio"/> 26"				<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$262.00	
<b>Back:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$169.11	<b>Back Extension:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$125.05
13" Height	<input type="radio"/> 22"	<input type="radio"/> 24"		\$191.73	12" Height	<input type="radio"/> 22"	<input type="radio"/> 24"		\$135.76
	<input type="radio"/> 26"			\$202.46		<input type="radio"/> 26"			\$147.67

<input type="radio"/> <b>Acute Care Configuration:</b>					
Includes: Full Vinyl Covers, Acute Care Short Shoulder Supports, ABS Lower Leg Support with Foot Support Cushions					
(Factory installed)(skip lower leg support cushion package section above when selecting this option)(16"-24" widths)					
<b>Common Accessories</b>					

<b>Neck Supports:</b>	<input type="radio"/> NR-CV	Neck Support Visco Foam Contoured Cushion	<input type="radio"/> Triptex	<input type="radio"/> Omit Standard	\$122.66
	<input type="radio"/> NS	Neck Support Soft Visco Cushion		Neck Support	\$122.66
<b>Positioning:</b>	<input type="radio"/> SR	Width Reduction Cushions - Pair			\$116.54
	<input type="radio"/> SAS	Side Arm Supports - Pair			\$143.58
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports - Pair			\$547.09
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED) - Pair			\$674.66
<b>Foot Cushions:</b>	<input type="radio"/> FB-LONG	Lower Leg Support Footbox - Full Extended Cushion	<input type="radio"/> Omit Foot Support Cushions		\$398.66
<b>Pelvic Belts:</b>	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned			\$88.32
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned & Auto Buckle			\$88.32

<b>Thigh Belts:</b>	<input type="radio"/> T-BELT	Thigh Belt Cushioned			\$141.06
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned & Auto Buckle			\$147.19
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene			\$141.06
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene Auto Buckle			\$147.19
	<input type="radio"/> T-BELT-NEO-APP	Thigh Belt Cushioned Neoprene for APP & SPP Packages			\$141.06
	<input type="radio"/> T-BELT-NEO-A-APP	Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages			\$147.19

<b>Trays:</b>	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$361.86
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)		<input type="radio"/> Black Only	\$521.32
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$374.13
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover			\$184.01

<b>Stabilizers:</b>	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one)(pair)			\$104.26
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers			\$184.01

<b>Other:</b>	<input type="radio"/> IV-POLE*	IV Pole and Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount	\$116.54
	<input type="radio"/> O2-HOLDER*	Oxygen Holder and Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount	\$196.27

\*IV-POLE and O2-HOLDER not compatible with rear mounted mag wheels

<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each	\$49.06
<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font OR <input type="radio"/> Block Font	\$30.67
<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag	\$104.94
<input type="radio"/> CUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)	\$55.19

**Additional Options:**

**Delivery Instructions or Special Requests:**

**Total: \$**

**Accessorial Services:** Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

**any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

**Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**