

PO:	Date:	Customer #:
Customer Name:	Ship To Name:	
Address:	Address:	
Phone:	Phone:	
Contact:	Contact:	
Email:	Tag:	



Standard Features	
· Tubular Steel, Powder Coated Frame	· Drop In, Swing Away Lower Leg Supports
· 45 Degrees, Infinitely Adjustable Seat Tilt	· Composite Foot Support Plates
· Fixed Height Arm Supports	· Cushioned Side Panels
· Adjustable Seat Depth 16" - 20"	· Adjustable Backrest Angle
· Adjustable Seat Surface Height 14" - 19"	· 250 lbs. Maximum Occupant Weight
· Choice of Wheel, Caster and Fork Sizes	· Seat, back and head support cushions not included
· Push Lock Parking Brakes	· Angle Adjustable Aluminum Foot Plates WAMJ

Comfort Rehab Wheelchair (587) - MSRP \$2,595


ADP Code: WMA5B-025

Frame Color Options

Frame Color: Black Blue Bronze Nickel Red Orange Sand Yellow Pink Purple

Seat and Back Options	
Seat Width:	Seat Depth:
<input type="radio"/> 16" <input type="radio"/> 18" N/C	<input type="radio"/> 16" <input type="radio"/> 17" N/C
<input type="radio"/> 20" WAMA \$195.00	<input type="radio"/> 18" <input type="radio"/> 19" WAMB N/C
<input type="radio"/> 22" WAMA \$195.00	<input type="radio"/> 20" WAMB N/C
<input checked="" type="radio"/> Solid Seat Pan Factory Installed N/C	

Seat Surface Height:	<input type="radio"/> 14" <input type="radio"/> 14.5" <input type="radio"/> 15" <input type="radio"/> 15.5" N/C	Back Cane	<input type="radio"/> Ultra-Low Double Bend (16") Fixed Standard
<input type="radio"/> 16" <input type="radio"/> 16.5" <input type="radio"/> 17" <input type="radio"/> 17.5" N/C		Height:	<input type="radio"/> Ultra Low (16") Fold Down \$220
<input type="radio"/> 18" <input type="radio"/> 19" N/C			<input type="radio"/> Low (18") Fold Down \$220
			<input type="radio"/> Medium (21") Fold Down \$220
			<input type="radio"/> Tall (24") Fold Down \$220



Refer to Seat Height & Wheel Combinations Matrix for assistance. □

Front Mount Options	Wheel Options	Rear Mount Options
<input type="radio"/> 5" Fixed Caster Wheels - Pair Standard	<input type="radio"/> 20" Mag Wheels - Pair Standard	
<input type="radio"/> 6" Fixed Caster Wheels - Pair N/C	<input type="radio"/> 22" Mag Wheels - Pair N/C	
<input type="radio"/> 8" Fixed Caster Wheels - Pair N/C	<input type="radio"/> 24" Mag Wheels - Pair N/C	

Caster Forks: Required	Mag Wheel Position: Factory Set (Adjustable in the future)
<input type="radio"/> Small - Recommended for Small Mag Wheel (20" or 22") and Lower Seat Surface Height N/C	<input type="radio"/> A (Most Forward) - Smaller Turning Radius N/C
<input type="radio"/> Medium - Recommended for Large Mag Wheel (24") Higher Seat Surface Heights N/C	<input type="radio"/> B N/C
	<input type="radio"/> C N/C
	<input type="radio"/> D* (Most Rearward) - Larger Turning Radius N/C
	*Do not recommended with 20" Mag Wheels

Refer to Seat Height & Wheel Combinations Matrix for assistance. □

Parking Brake Extension Handles - Pair N/C

Anti-Tip Devices - Required Safety Features - Factory Installed N/C

Arm Support Options

Height Adjustable, Removable Arm Supports N/C

Lower Leg Support Assembly and Lower Leg Support Options	
Lower Leg Support Assembly Style: <i>must choose one</i>	Lower Leg Support Accessories:
<input type="radio"/> Drop In Straight - Pair N/C	<input type="radio"/> Posterior Foot Supports - Pair \$20
<input type="radio"/> Drop In Flared - Pair N/C	<input type="radio"/> Lower Leg Support Panel (for calf support) \$60
<input type="radio"/> Pin and Saddle Straight - Pair N/C	<input type="radio"/> Lower Leg Support Foot Board \$110
<input type="radio"/> Pin and Saddle Flared - Pair N/C	<input type="radio"/> F-SW-ELV Swing Away, Elevating Lower Leg Supports Black - Pair WAMK \$325
<input type="radio"/> Drop In IKON Style - Pair N/C	<input type="radio"/> F-SW-ELV-IKN Swing Away, Elevating Lower Leg Supports Silver - Pair WAMK \$405
<input type="radio"/> F-SW75ADJ-587V3 Swing Away 75° Lower Leg Supports - Pair N/C	<input type="radio"/> F-SW-ANGLE Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair WAMJ \$284
<input type="radio"/> F-SW60ADJ-587V3 Swing Away 60° Lower Leg Supports - Pair N/C	<input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support w/ Angle Adjustable Footplates Black - Pair WAMJ \$284
<input type="radio"/> F-SW-IKN Swing Away Footrests Lower Leg Supports Silver - Pair N/C	
<input type="radio"/> F-SW-IKN-ADJ Swing Away Lower Leg Supports Black w/ Angle Adjustable Flip-up Foot Supports N/C WAMJ	

Common Accessories

Pelvic Strap Pelvic Positioning Strap Grommet Mount Ctr Buckle Side Buckle \$52

Trays:

ABS Tray with Strap (16", 18" 20" SETND0005, 22" SETND0010 - Matches Corresponding Seat Width) Sewn Strap Grey Black \$304

Padded ABS Tray Cover (ABS Tray only) SETND0035 \$155

Clear Tray with Strap (16", 18", 20" - Matches Corresponding Seat Width) Sewn Strap \$314

Other:

IV Pole and Assembly Front Mounted OR Rear Mounted \$98

Oxygen Holder and Assembly WAMX \$165

Accessory Bag (Black) Each \$41

Accessory Bag Personalizing 2 Initials Brush Font OR Block Font \$26

O2-HOLDER-BAG Oxygen Holder, Soft Carry Bag \$99.00

Cup Holder \$46

Matrx® Seating Series
Matrx® Seat Cushions



For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth: 16" 18" 20" 22" (Only available w/ PS Cushion)

Seat Cushion Width: 16" 18" 20" 22"

<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown	\$797.15	<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown	\$1,025
<input type="radio"/> Multi - Recommended for moderate to high risk of skin breakdown	\$745	<input type="radio"/> Multi - Recommended for moderate to high risk of skin breakdown	\$995
<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown	\$667.68	<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown	\$681
<input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown	\$690.15	<input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown	
<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown	\$608.83	<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown	\$975

Libra Cushion Accessories:	<input type="radio"/> Libra Fluid Supplement	\$120	Vi, PSP & PS Cushion Accessories:	<input type="radio"/> Pelvic Obliquity Kit	\$85
	<input type="radio"/> Libra 1" Hemi Wedge	\$101			
	<input type="radio"/> Libra 1" Anterior Seat Wedge	\$202			
	<input type="radio"/> Libra Tapered Wedge Left	\$85	Cushion Rigidizer:		
	<input type="radio"/> Libra Tapered Wedge Right	\$85	All Cushions		
	<input type="radio"/> Libra Abductor	\$85			
				16" Depth	18" Depth
				16" Width	<input type="radio"/> \$115
				18" Width	<input type="radio"/> \$115
				20" Width	<input type="radio"/> \$115
				22" Width	N/A

Matrx® Backs

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height: 16" 18" (excludes PB & PB Deep) 20"

Back Width: 16" 18" 20" 22"

<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications	\$866	<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications	\$1,151
<input type="radio"/> Elite - Recommended for more active users	\$766	<input type="radio"/> Elite - Recommended for more active users	\$1,051
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs	\$918	<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs	\$1,095
<input type="radio"/> PB - Recommended for less active users seeking comfort	\$734		
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support	\$892		

Lateral Support Accessories:	<input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L	\$116	Privacy Flap Accessories:	<input type="radio"/> 16" Width Small Privacy Flap	\$53
	<input type="radio"/> Elite Back Swing away Lateral Hardware	\$243		<input type="radio"/> 18" Width Medium Privacy Flap	\$53
	<input type="radio"/> Swing away Lateral Pad (each) 4.25"Hx6.5"L	\$122		<input type="radio"/> 20" Width Large Privacy Flap	\$53
				<input type="radio"/> 22" Width Extra Large Privacy Flap	\$53

Matrx® Headrest

Head Support Mounting:	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support	\$285
Choose One	<input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support	\$275
Head Support Options:	<input type="radio"/> Elan Head Support	\$310
	<input type="radio"/> Elan 4 Point Head Support	\$448
	<input type="radio"/> Elan Headband with 4 Point Cushion	\$610

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorial services needed to deliver are requested when an order is submitted to Broda.

Lift Gate Service: Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

Residential Delivery Service: Delivery to a residential location

Notification Before Delivery: Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

Inside Delivery Service: Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

Limited Access Service: Delivery locations that are not open to the walk-in public, during normal business hours, Mon-Fri, 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.

587V3: SEAT HEIGHT & WHEEL COMBINATIONS								
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB
14"	20	5	SMALL	1	1	TOP	3	2
15"	20	5	SMALL	2	2	TOP	2	1
	22	5	SMALL	1	2	TOP	2	1
15.5"	20	6	MEDIUM	2	1	TOP	2	1
	22	6	MEDIUM	1	1	TOP	2	1
	24	6	SMALL	1	2	TOP	1	1
	24	5	MEDIUM	1	2	TOP	1	1
16"	24	6	MEDIUM	1	1	TOP	1	1
	20	5	SMALL	3	3	TOP	1	1
	20	5	MEDIUM	3	3	TOP	1	1
	20	6	MEDIUM	3	2	TOP	1	1
	22	5	SMALL	2	3	TOP	1	1
	22	6	SMALL	2	2	TOP	1	1
	22	5	MEDIUM	2	3	TOP	1	1
16.5"	22	6	MEDIUM	2	2	TOP	1	1
	24	5	SMALL	1	3	TOP	1	1
	20	6	SMALL	3	3	TOP	1	1
	24	6	SMALL	2	3	BOTTOM	3	1
17"	24	5	MEDIUM	2	4	BOTTOM	3	1
	24	6	MEDIUM	2	3	BOTTOM	3	1
	20	5	MEDIUM	4	5	BOTTOM	3	1
	20	6	MEDIUM	4	4	BOTTOM	3	1
17.5"	22	5	MEDIUM	3	5	BOTTOM	3	1
	22	6	MEDIUM	3	4	BOTTOM	3	1
	20	8	SMALL	4	3	BOTTOM	3	1
	22	8	SMALL	3	3	BOTTOM	3	1
18"	24	8	SMALL	3	3	BOTTOM	2	1
	24	6	MEDIUM	3	5	BOTTOM	2	1
	24	8	MEDIUM	3	3	BOTTOM	2	1
	20	8	MEDIUM	5	4	BOTTOM	2	1
19"	22	8	MEDIUM	4	4	BOTTOM	2	1
	24	8	MEDIUM	4	5	BOTTOM	1	1

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.

Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG- 404 587V3 Rev. 1