

PO:	Date:	Customer #:
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
· Tubular Steel, Powder Coated Frame	· Drop In, Swing Away Lower Leg Supports
· 45 Degrees, Infinitely Adjustable Seat Tilt	· Composite Foot Support Plates
· Fixed Height Arm Supports	· Cushioned Side Panels
· Adjustable Seat Depth 16" - 20"	· Adjustable Backrest Angle
· Adjustable Seat Surface Height 14" - 19"	· 250 lbs. Maximum Occupant Weight
· Choice of Wheel, Caster and Fork Sizes	· Seat, back and head support cushions not included
· Push Lock Parking Brakes	

Comfort Rehab Wheelchair (587) - MSRP \$2,888.24


Frame Color Options

Frame Color: Black Blue Bronze Nickel Red Orange Sand Yellow Pink Purple

Seat and Back Options

Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	N/C	Seat Depth:	<input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20"	N/C
		\$206.70			N/C
		\$206.70			N/C
	<input checked="" type="radio"/> Solid Seat Pan	Factory Installed			N/C

Seat Surface Height:	<input type="radio"/> 14" <input type="radio"/> 14.5" <input type="radio"/> 15" <input type="radio"/> 15.5" <input type="radio"/> 16" <input type="radio"/> 16.5" <input type="radio"/> 17" <input type="radio"/> 17.5" <input type="radio"/> 18" <input type="radio"/> 19"	N/C	Back Cane Height:	<input type="radio"/> Ultra-Low Double Bend (16") Fixed <input type="radio"/> Ultra Low (16") Fold Down <input type="radio"/> Low (18") Fold Down <input type="radio"/> Medium (21") Fold Down <input type="radio"/> Tall (24") Fold Down	Standard \$233.20 \$233.20 \$233.20 \$233.20
-----------------------------	---	-----	--------------------------	---	--



Refer to Seat Height & Wheel Combinations Matrix for assistance. □

Back Cane Angle: 5° 10° 15° 20°

Front Mount Options

Wheel Options

Rear Mount Options

<input type="radio"/> 5" Fixed Caster Wheels - Pair	Standard	<input type="radio"/> 20" Mag Wheels - Pair	Standard
<input type="radio"/> 6" Fixed Caster Wheels - Pair	N/C	<input type="radio"/> 22" Mag Wheels - Pair	N/C
<input type="radio"/> 8" Fixed Caster Wheels - Pair	N/C	<input type="radio"/> 24" Mag Wheels - Pair	N/C

Caster Forks: Required		Mag Wheel Position: Factory Set (Adjustable in the future)	
<input type="radio"/> Small - Recommended for Small Mag Wheel (20" or 22") and Lower Seat Surface Heights	N/C	<input type="radio"/> A (Most Forward) - Smaller Turning Radius	N/C
<input type="radio"/> Medium - Recommended for Large Mag Wheel (24") Higher Seat Surface Heights	N/C	<input type="radio"/> B	N/C
		<input type="radio"/> C	N/C
		<input type="radio"/> D* (Most Rearward) - Larger Turning Radius	N/C
		*Do not recommended with 20" Mag Wheels	

Refer to Seat Height & Wheel Combinations Matrix for assistance. □

Parking Brake Extension Handles - Pair N/C Anti-Tip Devices - Required Safety Features - Factory Installed N/C

Arm Support Options

Height Adjustable, Removable Arm Supports N/C

Lower Leg Support Assembly and Lower Leg Support Options

Lower Leg Support Assembly Style: must choose one		Lower Leg Support Accessories:	
<input type="radio"/> Drop In Straight - Pair	<input type="radio"/> Drop In Flared - Pair	<input type="radio"/> Posterior Foot Supports - Pair	\$21.20
<input type="radio"/> Pin and Saddle Straight - Pair	<input type="radio"/> Pin and Saddle Flared - Pair	<input type="radio"/> Lower Leg Support Panel (for calf support)	\$63.60
<input type="radio"/> Drop In IKON Style - Pair		<input type="radio"/> Lower Leg Support Foot Board	\$116.60
<input type="radio"/> F-SW75ADJ-587V3 Swing Away 75° Lower Leg Supports - Pair	N/C	<input type="radio"/> F-SW-ELV Swing Away, Elevating Lower Leg Supports Black - Pair	\$344.50
<input type="radio"/> F-SW60ADJ-587V3 Swing Away 60° Lower Leg Supports - Pair	N/C	<input type="radio"/> F-SW-ELV-IKN Swing Away, Elevating Lower Leg Supports Silver - Pair	\$429.30
<input type="radio"/> F-SW-IKN	N/C	<input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support w/ Angle Adjustable Footplates Black - Pair	\$429.30
<input type="radio"/> F-SW-IKN-ADJ Swing Away Lower Leg Supports Black w/ Angle Adjustable Flip-up Foot Supports	N/C		
<input type="radio"/> 2" Extended Footrest Tubes - Only compatible on the F-SW-IKN and F-SW-IKN-ADJ	\$71.45		

Common Accessories

Pelvic Strap	<input type="radio"/> Pelvic Positioning Strap Grommet Mount	<input type="radio"/> Ctr Buckle	<input type="radio"/> Side Buckle	\$55.12
Trays:	<input type="radio"/> ABS Tray with Strap (16", 18" 20" 22" - Matches Corresponding Seat Width)	<input type="radio"/> Sewn Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$361.86
	<input type="radio"/> Padded ABS Tray Cover (ABS Tray only)			\$184.01
	<input type="radio"/> Clear Tray with Strap (16", 18", 20" - Matches Corresponding Seat Width)	<input type="radio"/> Sewn Strap		\$374.13
Other:	<input type="radio"/> IV Pole and Assembly	<input type="radio"/> Front Mounted	OR <input type="radio"/> Rear Mounted	\$116.54
	<input type="radio"/> Oxygen Holder and Assembly			\$196.27
	<input type="radio"/> Accessory Bag (Black) Each			\$49.06
	<input type="radio"/> Accessory Bag Personalizing 2 Initials	<input type="radio"/> Brush Font	OR <input type="radio"/> Block Font	\$30.67
	<input type="radio"/> O2-HOLDER-BAG Oxygen Holder, Soft Carry Bag			\$104.94
	<input type="radio"/> Cup Holder			\$55.19

Matrx® Seating Series

Please note, when ordering the Seat cushion only please contact Broda for additional shipping costs



Matrx® Backs - All Back Support Options Ship with Mounting Hardware Pre-Configured on Back Canes

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height: <input type="radio"/> 16" <input type="radio"/> 18" (excludes PB & PB Deep) <input type="radio"/> 20"	
Back Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" Back Width: <input type="radio"/> 22"	
<input type="radio"/> E2 - Recommended for users seeking optimal immersion and airflow \$965.00	<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$1,151.00
<input type="radio"/> E2 Deep - Recommended for users seeking optimal immersion w/ contoured support needs \$1,125.00	<input type="radio"/> Elite - Recommended for more active users \$1,124.57
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$866.00	<input type="radio"/> Elite Deep - Recommended for active users with contoured support \$1,171.65
<input type="radio"/> Elite - Recommended for more active users \$766.00	Matrx Back <input type="radio"/> Air Force N/C
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs \$918.00	Military <input type="radio"/> Army N/C
<input type="radio"/> PB - Recommended for less active users seeking comfort \$734.00	Service <input type="radio"/> Marine N/C
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support \$892.00	Decals: <input type="radio"/> Navy N/C
E2 Back Accessories:	Privacy Flap Accessories:
<input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair) \$385.00	<input type="radio"/> 16" Width Small Privacy Flap \$56.71
<input type="radio"/> E2 Chest Harness Interface \$85.00	<input type="radio"/> 18" Width Medium Privacy Flap \$56.71
<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support \$75.00	<input type="radio"/> 20" Width Large Privacy Flap \$56.71
<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32"L \$295.00	<input type="radio"/> 22" Width Extra Large Privacy Flap \$56.71
<input type="radio"/> Quad Release Tab for PostureFlex \$35.00	Lateral Support Accessories:
<input type="radio"/> E2 Infection Control Slipcover \$95.00	<input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L \$124.12
<input type="radio"/> E2 Deep Infection Control Slipcover \$95.00	<input type="radio"/> Elite Back Swing Away Lateral Hardware \$260.01
<input type="radio"/> Infection Control Slipcover for Matrx Cushions \$95.00	<input type="radio"/> Swing Away Lateral Pad (Each) 4.25"Hx6.5"L \$130.54

Matrx® Seat Cushions

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" (Only available w/Multi & PS Cushions)																					
Seat Cushion Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" Seat Cushion Width: <input type="radio"/> 22"																					
<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown \$797.15	<input type="radio"/> Libra - Recommended for highest risk or existing skin breakdown \$1,096.75																				
<input type="radio"/> Multi - Recommended for moderate to high risk of skin breakdown \$745.00	<input type="radio"/> Multi - Recommended for moderate to high risk of skin breakdown \$995.00																				
<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$667.68	<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$728.67																				
<input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown \$690.15	<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$1,043.25																				
<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$608.83	Multi, Vi, PSP & PS Cushion Accessories: <input type="radio"/> Pelvic Obliquity Kit \$90.95																				
Libra Cushion Accessories:	Cushion Rigidizer:																				
<input type="radio"/> Libra Fluid Supplement \$128.40	<table border="1"> <thead> <tr> <th></th> <th>16" Depth</th> <th>18" Depth</th> <th>20" / 22" Depth</th> </tr> </thead> <tbody> <tr> <td>All Cushions 16" Width</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> <tr> <td>18" Width</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> <tr> <td>20" Width</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> <tr> <td>22" Width</td> <td>N/A</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> </tbody> </table>		16" Depth	18" Depth	20" / 22" Depth	All Cushions 16" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	18" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	20" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	22" Width	N/A	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05
		16" Depth	18" Depth	20" / 22" Depth																	
All Cushions 16" Width		<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																	
18" Width		<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																	
20" Width		<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																	
22" Width	N/A	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																		
<input type="radio"/> Libra 1" Hemi Wedge \$108.07																					
<input type="radio"/> Libra 1" Anterior Seat Wedge \$216.14																					
<input type="radio"/> Libra Tapered Wedge Left \$90.95																					
<input type="radio"/> Libra Tapered Wedge Right \$90.95																					
<input type="radio"/> Libra Abductor \$90.95																					

Matrx® Headrest

Head Support Mounting:	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support \$285
Choose One	<input type="radio"/> Elan Flipp Head Support Hardware - Required with Selection of Head Support \$715
	<input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support \$275
	<input type="radio"/> Loxx Flipp Head Support Hardware - Required with Selection of Head Support \$705
Support Options:	<input type="radio"/> Head Support Cushion <input type="radio"/> 10" <input type="radio"/> 14" \$331.70
	<input type="radio"/> Hug Head Support Cushion <input type="radio"/> 10" <input type="radio"/> 14" \$331.70
	<input type="radio"/> Head Support 4 Point Cushion <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10" \$479.36
Additional Options:	<input type="radio"/> Headband with 4 Point Cushion <input type="radio"/> Small <input type="radio"/> Large \$652.70

Delivery Instructions or Special Requests:

Total: \$

Accessory Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. To ensure a seamless delivery of each Broda order,

any accessory services needed to deliver are requested when an order is submitted to Broda.

Lift Gate Service: Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

Residential Delivery Service: Delivery to a residential location

Notification Before Delivery: Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

Inside Delivery Service: Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

Limited Access Service: Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.

Knee-to-Heel Range (Comfort Rehab ONLY)			
Front Rigging Model	Description	Minimum	Maximum
Range Adjusts in 3/4" increments			
F-SW 75-587V3	75° Drop-in, black	14.25"	17.25"
F-SW 75ADJ-587V3	75° Drop-in w/ angle adjustable footplates, black	14.25"	17.25"
F-SW 60ADJ-587V3	60° Drop-in w/ angle adjustable footplates, black	15"	18"
F-SW-ELV-587V3	Drop-in, elevating, black	15"	18"
F-SW-ELV-ANGLE-587V3	Drop-in, elevating w/ angle adjustable footplates, black	15"	18"

587V3: SEAT HEIGHT & WHEEL COMBINATIONS								
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB
14"	20	5	SMALL	1	1	TOP	3	2
	20	5	SMALL	2	2	TOP	2	1
15"	22	5	SMALL	1	2	TOP	2	1
	20	6	MEDIUM	2	1	TOP	2	1
15.5"	22	6	MEDIUM	1	1	TOP	2	1
	24	6	SMALL	1	2	TOP	1	1
	24	5	MEDIUM	1	2	TOP	1	1
	24	6	MEDIUM	1	1	TOP	1	1
16"	20	5	SMALL	3	3	TOP	1	1
	20	5	MEDIUM	3	3	TOP	1	1
	20	6	MEDIUM	3	2	TOP	1	1
	22	5	SMALL	2	3	TOP	1	1
	22	6	SMALL	2	2	TOP	1	1
	22	5	MEDIUM	2	3	TOP	1	1
	22	6	MEDIUM	2	2	TOP	1	1
16.5"	24	5	SMALL	1	3	TOP	1	1
	20	6	SMALL	3	3	TOP	1	1
	24	6	SMALL	2	3	BOTTOM	3	1
	24	5	MEDIUM	2	4	BOTTOM	3	1
17"	24	6	MEDIUM	2	3	BOTTOM	3	1
	20	5	MEDIUM	4	5	BOTTOM	3	1
	20	6	MEDIUM	4	4	BOTTOM	3	1
	22	5	MEDIUM	3	5	BOTTOM	3	1
17.5"	22	6	MEDIUM	3	4	BOTTOM	3	1
	20	8	SMALL	4	3	BOTTOM	3	1
	22	8	SMALL	3	3	BOTTOM	3	1
	24	8	SMALL	3	3	BOTTOM	2	1
18"	24	6	MEDIUM	3	5	BOTTOM	2	1
	24	8	MEDIUM	3	3	BOTTOM	2	1
	20	8	MEDIUM	5	4	BOTTOM	2	1
	22	8	MEDIUM	4	4	BOTTOM	2	1
19"	24	8	MEDIUM	4	5	BOTTOM	1	1

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.

Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG-404 587V3 Rev. 1