

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



(shown with optional upgrades)

Standard Features	
• 250lb Maximum Occupant Weight (rocking)	• Rocking Mechanism
• 350lb Maximum Occupant Weight (non-rocking)	• Seat Mount Lower Leg Support Assembly
• Tubular Steel, Powder Coated Frame	• 2 x 20" Front Mag wheels with Front Anti-Tippers
• Adjustable Seat Surface Height	• 2 x 5" Total Lock Rear Caster Wheels
• Independent Adjustable Arm Support Height	
• Up to 40° Infinitely Adjustable Seat Tilt (non-rocking)	
• Up to 33° Infinitely Adjustable Seat Tilt (rocking)	

Encore Pedal Wheelchair® - Model: 48V4/48RV4 MSRP \$3,890.89

<input type="radio"/> WC19 Transit Tie Down Vehicle Package* - Factory Installed Only	\$416.82
*Available up to 24" seat width unless selected with the Heavy Duty package.	
<input type="radio"/> Heavy Duty Package** - Factory Installed Only, Increases Maximum Occupant Weight to 500 lbs.	\$550.00
**Available in seat widths 22"-28" only. Non-rocking only. Aluminum front casters not available.	
Sub-mount seat hanger style only with pin & saddle, angle-adjustable footplates (F-SW-ANGLE) only.	
If combined with WC19 package, only mag wheels are available (20", 22" front or 14" rear with twin front casters)	

Seat and Back Support Options

Rocking Mechanism: <input type="radio"/> Eliminate Rocking Function (\$269.14)	Seat Depth: <input type="radio"/> 17" Standard <input type="radio"/> 18.5" N/C <input type="radio"/> 20" \$116.54
Seat Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C <input type="radio"/> 22" \$484.53	Back Support Height: <input type="radio"/> 32" Standard <input type="radio"/> 33.5" \$116.54 <input type="radio"/> 35" \$116.54
Non-Rocking Only: <input type="radio"/> 24" <input type="radio"/> 26" \$1,097.84 <input type="radio"/> 28" \$1,097.84	Short Back Support Height: <input type="radio"/> 22" (does not include and Assembly) (\$490.65) N/C Height: Neck Support Option - see page 2
Seat Surface Height Rocking: <input type="radio"/> 15" Standard <input type="radio"/> 14" <input type="radio"/> 16" <input type="radio"/> 17" N/C <input type="radio"/> 18" <input type="radio"/> 19" N/C	
Seat Surface Height No Rock: <input type="radio"/> 12.5" <input type="radio"/> 13.5" <input type="radio"/> 14.5" N/C <input type="radio"/> 15.5" <input type="radio"/> 16.5" <input type="radio"/> 17.5" N/C	

Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair Standard	
<input type="radio"/> 22" Mag Wheels - Pair N/C	<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels* \$233.06 <i>*(available with front mags only)</i>	
<input type="radio"/> 5" Total Lock Caster Wheels - Pair	<input type="radio"/> 14" Mag Wheels - Pair N/C	
<input type="radio"/> 5" Total Lock Caster Wheels - Pair N/C	<input type="radio"/> 5" Fixed Non-Locking Caster Wheels - Pair (\$644.28)	
<input type="radio"/> Parking Brake Extension Handles (pair) \$78.50	<input type="radio"/> Front Mount Wheel Locks - Pair N/C	

Arm Support Options			
Arm Support	<input type="radio"/> Arm Support - Pair	Standard	Arm Support
Cushion	<input type="radio"/> Arm Support Trough - Each	\$159.46	Style:
Options:	<input type="radio"/> Cushioned Flip Arm Support:	\$142.05	<input type="radio"/> Height Adjustable Arm Supports
	<input type="radio"/> Cushioned Arm Support* -Pr	\$92.00	<input type="radio"/> Height and Depth Adjustable -Pair
	<i>*Not compatible with a tray</i>		<input checked="" type="radio"/> Seat Mounted Arm Support - Pair
			Mount:

Seating and Frame Color Options			
Comfort Tension	<input checked="" type="radio"/> Black	Standard	Frame Color:
Seating Color:		N/C	
			<input type="radio"/> Tungsten Sand
			<input type="radio"/> Bronze
			<input type="radio"/> Deep Orange
			<input type="radio"/> Carmine Red
			<input type="radio"/> Ocean Blue
			<input type="radio"/> Passion Pink
			<input type="radio"/> Satin Nickel
			<input type="radio"/> Black
			<input type="radio"/> Royal Purple

Upholstery Style Options (Select one)			
Accent Color Vinyl:	<input type="radio"/> Interlace Indigo***	<input type="radio"/> Interlace Espresso***	<input type="radio"/> Interlace Graphite***
	<input type="radio"/> Modena Black		
Vinyl:	<input type="radio"/> Beeline Portobello	<input type="radio"/> Kaleidoscope Smokehouse	<input type="radio"/> Beeline Forge
			<input type="radio"/> Beeline Geyser

Upholstery Colors come with seat and back cushions in black vinyl. *Interlace Indigo, Espresso, and Graphite available with seat and back cushions in matching vinyl color

Cushions - Choose Option 1, 2, OR 3			
<input type="radio"/> Option 1 - Broda Standard Cushions			N/C
Includes: Neck Support, Adjustable Lateral Shoulder Supports, Lower Side Covers, Seat and Back Support Cushions with Covers			
Seat & Back Pad Cover	<input type="radio"/> Black Vinyl		Standard
Options:	<input type="radio"/> Accent Color Vinyl	***Option For Interlace: Indigo, Espresso or Graphite only	N/C
	(Select Black Vinyl OR Accent Color Vinyl Above)		
	<input type="radio"/> Full Vinyl Wrap Seat and Back Support Cushions		\$116.54
	(Select Black Vinyl OR Accent Color Vinyl Above)		
	<input type="radio"/> Grey Terry Cover	<input type="radio"/> Startex Cover	N/C
Lateral Shoulder Support	<input type="radio"/> Lateral Shoulder Supports - Pair		N/C
Options:	<input type="radio"/> Long Lateral Shoulder Supports - Pair		\$159.46

Optional Substitutions Available Between APP & SPP Components

<input type="radio"/> Option 2 - Broda Additional Positioning Padding (APP), 24"-28" seat widths for non-rocking only			
Broda APP Package Includes: Seat, Back, Head Support Cushions, and Lateral Shoulder Supports			
on short back), plus Adjustable, Removable Lateral Shoulder Supports and Lateral Upper Leg Support Cushions			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21" \$462.44
	<input type="radio"/> 22" <input type="radio"/> 24"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21" \$706.54
	<input type="radio"/> 26" <input type="radio"/> 28"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21" \$836.56

<input type="radio"/> Option 3 - Broda Slim Posture Padding (SPP), 24"-26" seat widths for non-rocking only			
Broda SPP Cushion Options: Seat, Back Support, and Back Support Extension Cushions			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20" \$211.98
	<input type="radio"/> 22" <input type="radio"/> 24"		<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20" \$233.42
	<input type="radio"/> 26"		<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20" \$262.00
Lower Back	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" \$169.11	Upper Back	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" \$125.05
Support:	<input type="radio"/> 22" <input type="radio"/> 24" \$191.73	Support:	<input type="radio"/> 22" <input type="radio"/> 24" \$135.76
13" Height	<input type="radio"/> 26" \$202.46	12" Height	<input type="radio"/> 26" \$147.67

Matrx® Head Support (Available with Short Back Option Only)

Head Support Mounting:	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support	\$285.00
Choose One	<input type="radio"/> Elan Flipp Head Support Hardware - Required with Selection of Head Support	\$715.00
	<input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support	\$275.00
	<input type="radio"/> Loxx Flipp Head Support Hardware - Required with Selection of Head Support	\$705.00
Support Options:	<input type="radio"/> Elan Head Support Cushion <input type="radio"/> 10" <input type="radio"/> 14"	\$331.70
	<input type="radio"/> Hug Head Support 4 Point Cushion <input type="radio"/> 10" <input type="radio"/> 14"	\$331.70
	<input type="radio"/> Elan Head Support 4 Point Cushion <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10"	\$479.36
	<input type="radio"/> Elan Head Support with 4 Point Cushion <input type="radio"/> Small <input type="radio"/> Large	\$652.70

Drop In Style Must Choose One Lower Leg Support Assembly Style Pin and Saddle Style

<input type="radio"/> Lower Leg Support Assembly (Drop-in)- Pair N/C	<input type="radio"/> Lower Leg Support Assembly (Pin & Saddle) - Pair*N/C
<i>*Pin & Saddle Lower Leg Support Assembly not available with 16" SW chair due to overlapping foot supports</i>	

Must Choose One Seat Mount Lower Leg Support Assembly Location


<input type="radio"/> Seat Frame Mount* N/C	<input type="radio"/> Sub-Seat Frame Mount** N/C
<i>*Not available on rocking version with 17" seat depth</i>	<i>**Improves the ease of the rocking mechanism</i>

**Not available on non-rocking version when choosing 12.5" seat surface height with pin & saddle style lower leg support assembly*

Seat Frame Mount



Sub-Seat Frame Mount



Drop In Style Lower Leg Support Options Pin and Saddle Style

<input type="radio"/> F-SW-IKN Swing Away Lower Leg Support Silver -Pair \$128.80	<input type="radio"/> F-SW Swing Away Lower Leg Supports Black - P \$196.27
<input type="radio"/> F-SW-ELV-IKN Angle Adjustable Lower Leg Supports \$533.58	<input type="radio"/> F-SW-ELV Elevating Swing Away Lower Leg Supports Black - Pair \$257.59
Silver Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair	
<input type="radio"/> F-SW-IKN-ADJ Swing Away Lower Leg Supports Black w/ Adjustable Flip-up Foot Supports - Pair \$239.37	<input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support w/ Angle Adjustable Footplates - Black \$300.52
<input type="radio"/> F-SW Swing Away Lower Leg Supports Black - Pair	<input type="radio"/> LEG-FOOT-ENCORE Lower Leg Support Foot Board \$275.99
	<i>*Available with F-SW & F-SW-ELV only</i>
<input type="radio"/> 2" Extended Footrest Tubes - Only compatible on the F-SW-IKN and F-SW-IKN-ADJ \$71.45	<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support ABS Cushions \$92.00

Encore Pedal with Seat Mounted Hangers	Compatible with 14.5" Seat Height and above
Encore Pedal with Sub-Seat Mounted Hangers	Compatible with 15.5" Seat Height and above
Encore Pedal with Rocker with Seat Mounted Hangers	Compatible with 15" Seat Height and above
Encore Pedal with Rocker with Sub-Seat Mounted Hangers	Compatible with 16" Seat Height and above

Common Accessories

Head	<input type="radio"/> NR-CV Head Support Visco Foam Contoured Cushion <input type="radio"/> Triptex <input type="radio"/> Omit	\$122.66
Supports:	<input type="radio"/> NS Neck Support Soft Visco Cushion <input type="radio"/> Standard Neck Support	\$122.66
Positioning	<input type="radio"/> SR Width Reduction Cushions - Pair	\$116.54
Supports:	<input type="radio"/> CC-LTS Comfort Company Lateral Supports - Pair	\$547.09
	<input type="radio"/> LTS Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)	\$674.66
Pelvic Belts:	<input type="radio"/> S-BELT Pelvic Belt <input type="radio"/> Center Buckle <input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-A Pelvic Belt Auto Buckle <input type="radio"/> CenterBuckle <input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-PD Pelvic Belt Cushioned	\$88.32
	<input type="radio"/> S-BELT-PD-A Pelvic Belt Cushioned and Auto Buckle	\$88.32

Thigh Belts:	<input type="radio"/> T-BELT	Thigh Belt Cushioned		\$141.06	
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned and Auto Buckle		\$147.19	
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene		\$141.06	
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene and Auto Buckle		\$147.19	
	<input type="radio"/> T-BELT-NEO-APP	Thigh Belt Cushioned Neoprene for APP & SPP Packages		\$141.06	
	<input type="radio"/> T-BELT-NEO-A-APP	Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages		\$147.19	
Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$361.86
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)		<input type="radio"/> Black Only	\$521.32
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$374.13
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes upper side pad)		<input type="radio"/> Grey <input type="radio"/> Black	\$496.79
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover			\$184.01
	Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one) (pair)		\$104.26
<input type="radio"/> STABIL-LAT		Lateral Frame Stabilizers		\$158.00	
Other:	<input type="radio"/> IV-POLE	IV Pole & Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount		\$116.54
	<input type="radio"/> O2-HOLDER	Oxygen Holder & Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount		\$196.27
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag	(Not recommended with 22" back)		\$104.94
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each			\$49.06
	<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font	<input type="radio"/> Block Font		\$30.67
	<input type="radio"/> CUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)			\$55.19

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorial services needed to deliver are requested when an order is submitted to Broda.

Lift Gate Service: Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

Residential Delivery Service: Delivery to a residential location

Notification Before Delivery: Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

Inside Delivery Service: Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

Limited Access Service: Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.