

<b>PO:</b> _____	<b>Date:</b> _____	<b>Customer #:</b> _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	




Standard Features	
<ul style="list-style-type: none"> <li>• Solid Seat Pan</li> <li>• Tubular Steel, Powder Coated Frame</li> <li>• Seat cushion, back &amp; head support not included</li> <li>• Independent Height Adjustable Arm Supports</li> <li>• Rocking Mechanism (<i>allows up to 33° of tilt</i>)</li> <li>• Fold Down and Angle Adjustable Back Canes</li> <li>• 33° Infinitely Adjustable Seat Tilt</li> </ul>	<ul style="list-style-type: none"> <li>• 250lb Maximum Occupant Weight</li> <li>• Adjustable Seat Surface Height</li> <li>• 20" Front Mount Mag Wheels</li> <li>• 2 x 5" Rear Total Lock Caster Wheels</li> <li>• Lower Leg Support Assembly</li> <li>• Removable Arm Supports</li> </ul>

**Encore Rehab Wheelchair® - Model: 500RW MSRP \$3,896.66**

**WC19 Transit Tie Down Vehicle Package\*** - Factory Installed Only \$416.82  
 \*Compatible only with F-SW-ANGLE footrest

Seat and Back Support Options			
<b>Seat Width:</b> <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	N/C N/C \$470.41	<b>Adjustable</b> <input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20"	N/C N/C \$113.13
<input checked="" type="radio"/> Solid Seat Pan		N/C	
<b>Back Cane</b> <input type="radio"/> Ultra-Low Double Bend (16") Fixed		Standard	
<b>Height:</b> <input type="radio"/> Low (18") Fold Down <input type="radio"/> Medium (21") Fold Down <input type="radio"/> Tall (24") Fold Down		N/C N/C N/C	



Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair	Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair
<input type="radio"/> 22" Mag Wheels - Pair	N/C	<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels - Pair
<input checked="" type="radio"/> Front Anti-Tip Devices Required Safety Fea	N/C	<input type="radio"/> Parking Brake Extension Handles - Pair
<input type="radio"/> Quick Release Axles - Pair	\$142.91	<input type="radio"/> Front Mount Wheel Locks - Pair
		<input type="radio"/> <b>Mag Wheel Position:</b> Factory Set (Adjustable in the future)
		<input type="radio"/> A (Most Forward) - Smaller Turning Radius
		<input type="radio"/> B
		<input type="radio"/> C
		<input type="radio"/> D* (Most Rearward) - Larger Turning Radius
		*Do not recommended with 20" Mag Wheels

Arm Support Options			
<b>Arm Support Options:</b> <input type="radio"/> Arm Supports - Pair	Standard	<b>Arm Support Style:</b> <input type="radio"/> Height Adjustable Arm Supports	Standard
<input type="radio"/> Arm Support Trough - Each	\$154.82	<input type="radio"/> *Height and Depth Adjustable Arm Supports	\$154.82
<input type="radio"/> Cushioned Flip Arm Supports	\$142.05	*Not compatible with 22" mag wheels in the 14", 15", & 16" seat heights	
<input type="radio"/> Cushioned Arm Supports* - Pair	\$89.32	<input checked="" type="radio"/> <b>Arm Support Mount:</b> Seat Mounted Arm Supports	Standard
*Not compatible with a tray			

Frame Color Options			
<b>Frame Color:</b> <input type="radio"/> Black <input type="radio"/> Orange <input type="radio"/> Carmine Red	<input type="radio"/> Tungsten Sand <input type="radio"/> Ocean Blue	<input type="radio"/> Bronze <input type="radio"/> Satin Nickel <input type="radio"/> Yellow Zinc <input type="radio"/> Passion Pink <input type="radio"/> Royal Purple	

Drop In Style	Must Choose Lower Leg Support Assembly Style	Pin and Saddle Style
<input type="radio"/> SM-HANGER-48RV4-RKR-IKN (pair)	N/C	<input type="radio"/> SM-HANGER-48RV4-RKR-PIN (pair) <i>Less \$54</i>
*Lower Leg Support Assembly not available with 16" SW chair due to overlapping footplates		

Drop In Style	Lower Leg Support Options	Pin and Saddle Style
<input type="radio"/> F-SW-IKl Swing Away Silver (pair)	\$125.05	<input type="radio"/> F-SW Swing Away Lower Leg Support Black - Pair
<input type="radio"/> F-SW-IKN ADJ Swing Away Lower Leg Su	\$239.37	<input type="radio"/> F-SW-ELV Elevating Swing Lower Leg Support
Black Pair w/ Angle Adjustable Flip-up Foot Supports		Black - Pair
<input type="radio"/> F-SW-ELV-IKN Elevating Swing Away Lower Leg S	\$518.04	<input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support
Silver w/ Angle Adjustable Flip-up Foot Supports		w/ Angle Adjustable Footplates Black - Pair
<input checked="" type="radio"/> Heel Loops - Pair (for 3 above drop-in style options)	\$37.10	<input type="radio"/> LEG-FOOT-ENCORE Lower Leg Support Foot Board
<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support	\$89.32	<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support Panel
*Slight ground interference may occur in the lowest seat height. Slightly tilt chair or elevate legrest as needed		

2" Extended Footrest Tubes - *Only compatible on the F-SW-IKN and F-SW-IKN-ADJ* \$71.45

Knee-to-Heel Range			
Front Rigging Model	Description	Minimum	Maximum
Range Adjusts in 1" increments			
F-SW-IKN	70° Drop-in, silver	13"	17" (19" with optional extender)
F-SW-IKN-ADJ	70° Drop-in w/ angle adjustable footplates, black	13"	17" (19" with optional extender)
Range Adjusts in 3/4" increments			
F-SW	70° Pin & Saddle, black	11"	14"
F-SW-ANGLE	70° Pin & Saddle w angle adjustable footplates, black	11"	14.75"
F-SW-ELV	Pin & Saddle, elevating, black	15.25"	18.25"
F-SW-ELV-ANGLE	Pin & Saddle, elevating w/ angle adjustable footplates, black	15.25"	18.25"
Infinite Range Adjustments			
F-SW-ELV-IKN	Drop-in, elevating w/ angle adjustable footplates, silver	13"	20"

"IKN" footrests are compatible with Encore Pedal, Encore Rehab, and Comfort Rehab Wheelchairs  
All other footrests are compatible with all Pedal and Rehab Wheelchairs

**Common Accessories**

<b>Pelvic Belts:</b>	<input type="radio"/> S-BELT Pelvic Belt	<input type="radio"/> Ctr Buck <input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-A Pelvic Belt Auto Buckle	<input type="radio"/> Ctr Buck <input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-587V3 Pelvic Belt Grommet Mount	<input type="radio"/> Ctr Buck <input type="radio"/> Side Buckle	\$61.34
<b>Trays:</b>	<input type="radio"/> TRAY-ABS Tray ABS with Strap	<input type="radio"/> Sewn Str <input type="radio"/> Grey <input type="radio"/> Black	\$361.86
	<input type="radio"/> TRAY-ABS-OVR Tray ABS with Strap (24" - 28" width)	<input type="radio"/> Black Only	\$521.32
	<input type="radio"/> TRAY-CLEAR Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap	\$374.13
	<input type="radio"/> TRAY-HALF Half Tray (side specific)(with upper side pad)	<input type="radio"/> Grey <input type="radio"/> Black	\$496.79
	<input type="radio"/> TRAY-COVER Tray ABS Cushion Cover		\$184.01
<b>Other:</b>	<input type="radio"/> IV-POLE IV Pole & Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount		\$116.54
	<input type="radio"/> O2-HOLDER O2 Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount		\$196.27
	<input type="radio"/> O2-HOLDER-BAG Oxygen Holder, Soft Carry Bag		\$104.94
	<input type="radio"/> STABIL-FR Rear Stabilizer Frame		\$104.26

**Common Cushion Upgrades**

<b>Positioning:</b>	<input type="radio"/> SR Width Reduction Cushions - Pair	\$116.54
---------------------	--	----------

**Matrx® Seating Series**

Please note, when ordering the Seat cushion only please contact Broda for additional shipping costs



**Matrx® Backs** - All Back Support Options Ship with Mounting Hardware Pre-Configured on Back Cane

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

<b>Back Height:</b>	<input type="radio"/> 16" <input type="radio"/> 18" (excludes PB & PB L) <input type="radio"/> 20"		
<b>Back Width:</b>	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"		
<input type="radio"/> E2 - Recommended for users seeking optimal immersion and airflow	\$965.00	<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications	\$1,151.00
<input type="radio"/> E2 Deep - Recommended for users seeking optimal immersion w/ contoured	\$1,125.00	<input type="radio"/> Elite - Recommended for more active users	\$1,124.57
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications	\$866.00	<input type="radio"/> Elite Deep - Recommended for active users with contoured support	\$1,171.65
<input type="radio"/> Elite - Recommended for more active users	\$766.00	<b>Matrx Back</b>	<input type="radio"/> Air Force N/C
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs	\$918.00	<b>Military</b>	<input type="radio"/> Army N/C
<input type="radio"/> PB - Recommended for less active users seeking comfort	\$734.00	<b>Service</b>	<input type="radio"/> Marine N/C
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support	\$892.00	<b>Decals:</b>	<input type="radio"/> Navy N/C
<b>E2 Back</b>	<input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair) \$385.00	<b>Privacy Flap</b>	<input type="radio"/> 16" Width Small Privacy Flap \$56.71
<b>Accessories</b>	<input type="radio"/> E2 Chest Harness Interface \$85.00		<input type="radio"/> 18" Width Medium Privacy Flap \$56.71
	<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support \$75.00		<input type="radio"/> 20" Width Large Privacy Flap \$56.71
	<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32" \$295.00		<input type="radio"/> 22" Width Extra Large Privacy Flap \$56.71
	<input type="radio"/> Quad Release Tab for PostureFlex \$35.00	<b>Lateral Support</b>	<input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L \$124.12
<b>Back</b>	<input type="radio"/> E2 Infection Control Slipcover \$95.00		<input type="radio"/> Elite Back Swing Away Lateral Hardware \$260.01
<b>Slipcovers:</b>	<input type="radio"/> E2 Deep Infection Control Slipcover \$95.00	<b>Accessories:</b>	<input type="radio"/> Swing Away Lateral Pad (Each) \$130.54
	<input type="radio"/> Infection Control Slipcover for Matrx Cushions \$95.00		4.25"Hx6.5"L

**Matrx® Seat Cushions**

*For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions*

<b>Seat Cushion Depth:</b> <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" (Only available w/ PS Cushion)																										
<b>Seat Cushion Width:</b> <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	<b>Seat Cushion Width:</b> <input type="radio"/> 22"																									
<input type="radio"/> <b>Libra</b> - Recommended for highest risk of or existing skin breakdown    \$797.15	<input type="radio"/> <b>Libra</b> - Recommended for highest risk or existing skin breakdown    \$1,096.75																									
<input type="radio"/> <b>Multi</b> - Recommended for moderate to high risk of skin breakdown    \$745.00	<input type="radio"/> <b>Multi</b> - Recommended for moderate to high risk of skin breakdown    \$995.00																									
<input type="radio"/> <b>Vi</b> - Recommended for moderate to high risk of skin breakdown    \$667.68	<input type="radio"/> <b>Vi</b> - Recommended for moderate to high risk of skin breakdown    \$728.67																									
<input type="radio"/> <b>PSP</b> - Recommended for moderate to high risk of skin breakdown    \$690.15																										
<input type="radio"/> <b>PS</b> - Recommended for low to moderate risk of skin breakdown    \$608.83	<input type="radio"/> <b>PS</b> - Recommended for low to moderate risk of skin breakdown    \$1,043.25																									
<b>Libra Cushion Accessories:</b>	<b>Vi, Multi, PSP &amp; PS Cushion Accessories:</b> <input type="radio"/> Pelvic Obliquity Kit    \$90.95																									
<input type="radio"/> Libra Fluid Supplement    \$128.40																										
<input type="radio"/> Libra 1" Hemi Wedge    \$108.07																										
<input type="radio"/> Libra 1" Anterior Seat Wedge    \$216.14																										
<input type="radio"/> Libra Tapered Wedge Left    \$90.95																										
<input type="radio"/> Libra Tapered Wedge Right    \$90.95																										
<input type="radio"/> Libra Abductor    \$90.95																										
	<b>Cushion Rigidizer:</b>																									
	<table border="1"> <thead> <tr> <th>All Cushions</th> <th>Column 1</th> <th>16" Depth</th> <th>18" Depth</th> <th>20" / 22" Depth</th> </tr> </thead> <tbody> <tr> <td>16" Width</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> <tr> <td>18" Width</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> <tr> <td>20" Width</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> <tr> <td>22" Width</td> <td>N/A</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> <td><input type="radio"/> \$123.05</td> </tr> </tbody> </table>	All Cushions	Column 1	16" Depth	18" Depth	20" / 22" Depth	16" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	18" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	20" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	22" Width	N/A	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05
All Cushions	Column 1	16" Depth	18" Depth	20" / 22" Depth																						
16" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																						
18" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																						
20" Width	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																						
22" Width	N/A	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05	<input type="radio"/> \$123.05																						

**Matrx® Headrest**

<b>Head Support Mounting:</b>	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support    \$285
<b>Choose One</b>	<input type="radio"/> Elan Flipp Head Support Hardware - Required with Selection of Head Support    \$715
	<input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support    \$275
	<input type="radio"/> Loxx Flipp Head Support Hardware - Required with Selection of Head Support    \$705
<b>Support Options:</b>	<input type="radio"/> Head Support Cushion <input type="radio"/> 10" <input type="radio"/> 14"    \$331.70
	<input type="radio"/> Hug Head Support Cushion <input type="radio"/> 10" <input type="radio"/> 14"    \$331.70
	<input type="radio"/> Head Support 4 Point Cushion <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10"    \$479.36
	<input type="radio"/> Headband with 4 Point Cushion <input type="radio"/> Small <input type="radio"/> Large    \$652.70

**Additional Options:** \_\_\_\_\_

**Delivery Instructions or Special Requests:** \_\_\_\_\_

Total: \$

**Accessorial Services:** Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

**any accessorial services needed to deliver are requested when an order is submitted to Broda.**

**Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

**Residential Delivery Service:** Delivery to a residential location

**Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

**Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

**Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

**Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**