

PO: _____ **Date:** _____ **Customer #:** _____

Customer Name: _____
Address: _____
Phone: _____
Contact: _____
Email: _____

Ship To Name: _____
Address: _____
Phone: _____
Contact: _____
Tag: _____



Standard Features	
• Comfort Tension Seating®	• 350lb Maximum Occupant Weight
• Up to 40° Infinitely Adjustable Seat Tilt	• Trendelenburg Capable
• Up to 90° Infinitely Adjustable Back Support Recline	• Angle Adjustable, Articulating Lower Leg Support
• 2 x 20" Rear Mag Wheels	• 5" Total Lock Front Caster Wheels
• Height Adjustable, Swing-Away Removable Arm Support	• Height and Width Adjustable Lateral Shoulder Supports

Synthesis Positioning Wheelchair - Model: V4 MSRP \$4,846.45

Seat and Back Support Options

Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C	Seat Depth:	<input type="radio"/> 17" Standard
	<input type="radio"/> 22" \$484.53		<input type="radio"/> 18.5" N/C
	<input type="radio"/> 24" <input type="radio"/> 26" \$1,097.84		<input type="radio"/> 20" \$116.54
			<input type="radio"/> 22" \$116.54
Seat Surface	<input type="radio"/> 19.5" \$208.53	Back Height:	<input type="radio"/> 32" N/C
Height:	<input type="radio"/> 20.5" Standard		<input type="radio"/> 33.5" \$116.54
	<input type="radio"/> 22" <input type="radio"/> 23.5" \$208.53		<input type="radio"/> 35" \$116.54

Front Options	Wheel Options	Rear Options
<input type="radio"/> 5" Total Lock Caster Wheels - Pair Standard	<input type="radio"/> 20" Mag Wheels - Pair Standard	<input type="radio"/> Omit Hand Rims N/C
<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels - Pair \$233.06	<input type="radio"/> 14" Mag Wheels - Pair N/C	<input type="radio"/> 6" Multi-Directional Locking Caster Wheels - Pair less \$385.85
<i>(available with rear mags only)</i>		
<input type="radio"/> 5" Locking Twin Caster Wheels - Pair \$233.06		

Middle Options

<input type="radio"/> 14" Mid Mags (comes with 5" Rear Caster Wheels)
<input type="radio"/> Parking Brake Extension Handles - Pair \$78.50

Lower Leg Support Options (Select One)

<input type="radio"/> Angle Adjustable, Articulating Split Lower Leg Support <i>(not available in 16" seat width)</i> Standard
<i>Choose Number of Straps (support length range)</i> <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Angle Adjustable, Articulating Single Lower Leg Support N/C
<i>Choose Number of Straps (support length range)</i> <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Angle Adjustable, Articulating ABS Lower Leg Support <i>(Includes calf and foot supports) (16" - 24" Widths)</i> \$398.66

Arm Support Options

Arm Support <input type="radio"/> Arm Supports - Pair Standard	Arm Support <input type="radio"/> Height Adjustable Arm Supports - Pair Standard
Cushion <input type="radio"/> Long Arm Supports - Pair \$49.06	Styles: <input type="radio"/> Height & Depth Adjustable Arm Sup \$159.46
Options: <input type="radio"/> Arm Support Trough - Each \$159.46	Arm Support <input checked="" type="radio"/> Seat Mounted Arm Supports Standard
<input type="radio"/> Cushioned Flip Arm Supports \$142.05	<input type="radio"/> Long Cushioned Arm Support* \$305.49
<input type="radio"/> Long Cushioned Flip Arm Support \$189.39	
<input type="radio"/> Cushioned Arm Supports* - Pair \$92.00	

*Not compatible with a tray

Seating Frame Color Options (Select one)

<input type="radio"/> Carmine Red	<input type="radio"/> Ocean Blue	<input type="radio"/> Midnight Black	N/C
<input type="radio"/> Satin Nickel	<input type="radio"/> Tungsten Sand	<input type="radio"/> Bronze	<input type="radio"/> Zinc Yellow N/C
<input type="radio"/> Passion Pink	<input type="radio"/> Royal Purple		N/C

Upholstery Color Options: (Select one)

- * Standard cushion configuration: Neck Support, Exterior Lateral Shoulder Supports, upper and lower side covers & foot cushions Accent Color. Seat, back and calf supports, and interior lateral shoulder supports come standard in Savoy Black Vinyl.
- ** Interlace: Indigo, Espresso & Graphite are available on seat, back & calf support cushions & most accessories.
- *** In honor of our Veterans' Service, Camo options are now available.



Standard Color Vinyl*: Shown with Beeline Geyser Upholstery & Tungsten Sand Frame Color

Accent Color Vinyl**: Shown with Interlace Indigo Upholstery & Black Frame Color *Chair Shown with ABS lower leg support

Upholstery Accent Color*: (Select one)

- Vinyl:**
- Interlace Indigo**
 - Interlace Graphite**
 - Kaleidoscope Smokehouse
 - Kaleidoscope Chocolate
 - Modena Black
 - Interlace Espresso**

- Vinyl:**
- Beeline Forge
 - Beeline Geyser
 - Beeline Fossil
 - Beeline Portobello

Nylon:***

- Retro Camo \$159.46

Cushions - Choose Option 1, 2, OR 3

Option 1 - Broda Standard Cushions

N/C

Standard Broda Cushions:

Includes: **In Accent Color:** Neck Support, Exterior Lateral Shoulder Supports, Upper & Lower Side Cover & Foot Support Cushion*.

In Black Vinyl: Seat & Back Support Cushions with Covers, Calf Support Cushion & Interior Lateral Shoulder Support Cushions

Seat, Back Support & Calf Support Cover Options:	<input type="radio"/> Black Vinyl	Standard
	<input type="radio"/> Accent Color*	N/C
	*** <i>(For Interlace: Indigo, Espresso or Graphite only)</i>	
	<input type="radio"/> Full vinyl Wrap Seat & Back Support Cushions <i>(Select Black Savoy Vinyl OR Accent Color Vinyl Above)</i>	\$116.54
	<input type="radio"/> Calf Support Cushion Color <i>(Select Black Savoy Vinyl OR Accent Color Vinyl Above)</i>	N/C
	<input type="radio"/> Grey Terry Cover <input type="radio"/> Startex Cover	N/C

Huntington's Special Padding (HSP) package:	<input type="radio"/> Includes Long Lateral Shoulder Supports, Thigh Belt & Long Lower Leg Support Footbox	\$631.72
	Thigh Belt Color Option: <input type="radio"/> Modena Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	Long Lower Leg Support Footbox Color: <input type="radio"/> Black Savoy Vinyl <input type="radio"/> Accent Color Vinyl	N/C

Optional Substitutions Available Between APP & SPP Components

Option 2 - Broda Additional Positioning Padding (APP) (not compatible with 35" back support height)

Broda APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Covers

Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$640.31
	<input type="radio"/> 22" <input type="radio"/> 24"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$770.33
	<input type="radio"/> 26"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$900.35

Option 3 - Broda Slim Posture Padding (SPP)

Broda SPP Padding Options: Seat, Back Support, and Back Support Extension Cushions

Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$211.98	
	<input type="radio"/> 22" <input type="radio"/> 24"		<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$233.42	
	<input type="radio"/> 26"		<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$262.00	
Back Support	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	\$169.11	Back Support	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	\$125.05
13" Height	<input type="radio"/> 22" <input type="radio"/> 24"	\$191.73	Extension:	<input type="radio"/> 22" <input type="radio"/> 24"	\$135.76
	<input type="radio"/> 26"	\$202.46	12" Height	<input type="radio"/> 26"	\$147.67

Acute Care Configuration:

Includes: Full Vinyl Covers, ABS Lower Leg Support with Foot Support Cushions (*Factory installed*)(skip lower leg support cushion \$552.00 package section above when selecting this option)(16"-24" widths)

Common Accessories					
Neck Supports:	<input type="radio"/> NR-CV	Head Support Visco Foam Contoured Cushion	<input type="radio"/> Triptex	<input type="radio"/> Omit Standard	\$122.66
	<input type="radio"/> NS	Neck Support Soft Visco Cushion		Neck Support	\$122.66
Positioning:	<input type="radio"/> SR	Width Reduction Cushions - Pair			\$116.54
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports - Pair			\$547.09
	<input type="radio"/> LTS	Lateral Torso Supports (<i>FACTORY INSTALLED RECOMMENDED</i>) - Pair			\$674.66
	<input type="radio"/> SB-LONG-V4	Long Lateral Shoulder Supports - Pair - Omits standard short			\$333.65
Foot Support Cushions:	<input type="radio"/> FB-LONG	Lower Leg Support Footbox - Full Extended Cushion		<input type="radio"/> Omit Foot Support Cushions	\$398.66
Pelvic Belts:	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned			\$88.32
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned & Auto Buckle			\$88.32
Thigh Belts:	<input type="radio"/> T-BELT	Thigh Belt Cushioned			\$141.06
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned & Auto Buckle			\$147.19
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene			\$141.06
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene Auto Buckle			\$147.19
	<input type="radio"/> T-BELT-NEO-APP	Thigh Belt Cushioned Neoprene for APP & SPP Packages			\$141.06
	<input type="radio"/> T-BELT-NEO-A-APP	Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages			\$147.19
Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$361.86
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 26" width)		<input type="radio"/> Black Only	\$521.32
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$374.13
	<input type="radio"/> TRAY-HALF	Half Tray (<i>side specific</i>)(includes upper side pad)			\$496.79
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover			\$184.01
Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (<i>select one</i>)(pair)			\$104.26
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers			\$184.01
Other:	<input type="radio"/> IV-POLE*	IV Pole and Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount		\$116.54
	<input type="radio"/> O2-HOLDER*	Oxygen Holder and Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount		\$196.27
	*IV-POLE and O2-HOLDER not compatible with rear mounted mag wheels				
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag			\$104.94
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each			\$49.06
	<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font OR <input type="radio"/> Block Font			\$30.67
	<input type="radio"/> CUP-HOLDER	Cup Holder			\$55.19

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

Lift Gate Service: Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

Residential Delivery Service: Delivery to a residential location

Notification Before Delivery: Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

Inside Delivery Service: Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

Limited Access Service: Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.