

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.



Standard Features	
• Comfort Tension Seating®	• 15° Adjustable Back Recline
• Tubular Steel, Powder Coated Frame	• 2 x 5" Total Lock Front Caster Wheels
• Infinitely Adjustable Seat Tilt to 40°	• 2 x 20" Mag Rear Wheels
• 350lb Maximum Occupant Weight	

Elite Rehab Wheelchair - Model: 550 SR (E1161) MSRP \$3,171.00

<input type="radio"/> WC-19 Transport Package: WC-19 Compliant Transit Tie Down Vehicle Transport Package - Factory Install Only	\$523.11
<input type="radio"/> Acute Care Configuration: Includes: Full Vinyl Covers, ACC Lateral Shoulder Supports, ABS Lower Leg Support with Foot Support Cushions (Factory Installed)(Skip lower leg support and cushion cover/wrap sections below when selecting this option, 16"-24" widths only)	\$545.37
<input type="radio"/> Huntington's Padding Package Includes: Long Lateral Shoulder Supports, Thigh Belt & Long Lower Foot Support Footbox	\$623.28

Seating Frame Color Options

Strap Color: <input type="radio"/> Grey <input type="radio"/> Black <input type="radio"/> Blue <input type="radio"/> N/C	Cushion <input type="radio"/> Burgundy <input type="radio"/> Milano Stitch Bark <input type="radio"/> N/C
Frame Color: <input type="radio"/> Grey <input type="radio"/> Black <input type="radio"/> N/C	Color: <input type="radio"/> Charcoal <input type="radio"/> Sapphire Blue <input type="radio"/> N/C
<input type="radio"/> Tungsten Sand <input type="radio"/> Ocean Blue	<input type="radio"/> Modena Black
<input type="radio"/> Carmine Red <input type="radio"/> Satin Nickel	
<input type="radio"/> Zinc Yellow <input type="radio"/> Deep Orange	
<input type="radio"/> Passion Pink <input type="radio"/> Royal Purple	
<input type="radio"/> Bronze	

Seat Options

Seat Width: <input type="radio"/> 16" or <input type="radio"/> 18" <input type="radio"/> N/C	Seat Depth: <input type="radio"/> 17" <input type="radio"/> Standard
<input type="radio"/> 20" (E2201) \$278.25	<input type="radio"/> 18.5" <input type="radio"/> N/C
<input type="radio"/> 22" (E2201) \$465.23	<input type="radio"/> 20" (E2203) \$278.25
<input type="radio"/> 24" (E2202) <input type="radio"/> 26" (E2202) \$1,052.90	<input type="radio"/> 22" (E2204) \$278.25
<input type="radio"/> 28" \$1,052.90	Seat Surface <input type="radio"/> 19" \$206.15
	Height: <input type="radio"/> 20.5" (ONLY option w/ WC-19 package) <input type="radio"/> Standard
	<input type="radio"/> 22" \$206.15
	<input type="radio"/> 23.5" \$206.15


Seat Option: <input type="radio"/> Comfort Tension Seating® (CTS) <input type="radio"/> Standard	<input type="radio"/> Solid Seat Pan (E2231) \$166.95
---	---

See CTS Seat Options below

See Matrix Seat Cushion Options below

Comfort Tension Seating® (CTS) ONLY: Seat Cushion Options

<input type="radio"/> Broda Seat Cushions Includes: <input type="radio"/> Vinyl Cover <input type="radio"/> Terry Cover (Grey) <input type="radio"/> Startex Cover <input type="radio"/> N/C	
<input type="radio"/> Full Vinyl Wrap Seat	\$57.32
<input type="radio"/> Broda Additional Positioning Padding (APP) Seat: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" E2601 Seat <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$188.91
Seat Cushion Width: <input type="radio"/> 22" <input type="radio"/> E2602 Depth: <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$213.44
<input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28" E2602 <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$252.69
<input type="radio"/> Broda Slim Posture Padding (SPP) Seat Cushion	
Seat Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" E2601 Seat Depth: <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$161.97
<input type="radio"/> 22" <input type="radio"/> 24" E2602 <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$178.64
<input type="radio"/> 26" E2602 <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	\$207.22

Solid Seat Pan: Matrx® ONLY: Seat Cushion Options			
For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions			
Seat Cushion Depth:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" (Only available w/Multi & PS Cushions)		
Seat Cushion Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24"(PS &Libra) <input type="radio"/> 26" (PS only)		
<input type="radio"/> Libra (E2624) -Recommend for highest risk of or existing skin breakdown	\$485	<input type="radio"/> Libra (E2625) -Recommend for highest risk of or existing skin breakdown	\$595
<input type="radio"/> Multi (E2607) -Recommend for moderate to high risk of skin breakdown	\$445	<input type="radio"/> Multi (E2608) -Recommend for moderate to high risk of skin breakdown	\$595
<input type="radio"/> Vi (E2607) -Recommend for moderate to high risk of skin breakdown	\$425	<input type="radio"/> Vi (E2608) -Recommend for moderate to high risk of skin breakdown	\$475
<input type="radio"/> PSP (E2607) -Recommend for moderate to high risk of skin breakdown	\$395		
<input type="radio"/> PS (E2605) - Recommend for low to moderate risk of skin breakdown	\$345	<input type="radio"/> PS (E2606) -Recommend for low to moderate risk of skin breakdown	\$545
Libra Cushion Accessories: <ul style="list-style-type: none"> <input type="radio"/> Libra Fluid Supplement \$75 <input type="radio"/> Libra 1" Hemi Wedge \$65 <input type="radio"/> Libra 1" Anterior Seat Wedge \$125 <input type="radio"/> Libra Tapered Wedge Left \$50 <input type="radio"/> Libra Tapered Wedge Right \$50 <input type="radio"/> Libra Abductor \$60 <input type="radio"/> Infection Control Slipcover for Matrx Cushions \$65 		Multi, Vi, PSP & PS Cushion Accessories: <ul style="list-style-type: none"> <input type="radio"/> Pelvic Obliquity Kit \$50 	
		Cushion Rigidizer:	16" D 18" D 20" / 22" D
		For All Cushions	16" Width <input type="radio"/> \$115 <input type="radio"/> \$115 <input type="radio"/> \$115
			18" Width <input type="radio"/> \$115 <input type="radio"/> \$115 <input type="radio"/> \$115
			20" Width <input type="radio"/> \$115 <input type="radio"/> \$115 <input type="radio"/> \$115
			22" Width N/A <input type="radio"/> \$115 <input type="radio"/> \$115

Back Options			
Back Height:	<input type="radio"/> 32" Standard <input type="radio"/> 33.5" <input type="radio"/> 35" \$111.30	Back Seating:	<input checked="" type="radio"/> Comfort Tension Seating® (CTS) Standard <input type="radio"/> Back Recline Full (E1226)(>15° - ≤ 80°) \$445.20
Back Recline:	<input type="radio"/> Back Recline Semi (0° - 15°) Standard		

Back Cushion Options			
<input checked="" type="radio"/> Broda Back Cushions Include: Upper & Lower Side Covers, Back Cushion with Cover, Foot Cushions			N/C
<input type="radio"/> Vinyl Back	<input type="radio"/> Terry Back (Grey)	<input type="radio"/> Startex Back	N/C
<input type="radio"/> Full Vinyl Wrap Back			\$57.32

Optional Substitutions Available Between APP & SPP Components					
<input type="radio"/> Broda Additional Positioning Padding (APP) Back Cushion		Back Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28"	E2611	\$179.10
				E2612	\$230.61
				E2612	\$257.60
<input type="radio"/> Broda Slim Posture Padding (SPP) Back Cushion		Back:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26"	E2611	\$107.19
13" Height		Back Extension	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26"	E2612	\$119.90
				E2612	\$131.10

Lateral Mount: <input type="radio"/> Adjustable Lateral Shoulder Support (E1028) (pair) (Upper lateral padding not included with lateral mount - select options below)		\$445.20
Lateral Shoulder Support Options: (lateral mount required)		
<input type="radio"/> Standard Short (E0956)(Pair)	\$148	<input type="radio"/> Rounded Short (E0956)(Pair) <input type="radio"/> ACC Short(E0956)(Pair) \$164.72
<input type="radio"/> Standard Long (E0956)(Pair)		<input type="radio"/> Rounded Long (E0956)(Pair) <input type="radio"/> ACC Long (E0956)(Pair) \$320.54
<input type="radio"/> APP V4 SB-SH	APP V4 Lateral Shoulder Supports - Pair \$122.66	
<input type="radio"/> Wing Covers Short(Pair)	\$199.23	
<input type="radio"/> Wing Covers Long (Pair)	(Wing covers can be used with standard short or standard long bolsters only) \$261.56	
Neck Support Option:	<input type="radio"/> NR (E0955) Neck Support	\$178.08
	<input type="radio"/> NR-CV (E0955) Neck Support Visco Foam Contoured Cushion <input type="radio"/> Triptex	\$121.32
	<input type="radio"/> NS (E0955) Neck Support Soft Visco Cushion	\$121.32
	<input type="radio"/> Additional Positionint Padding (APP) Head Support	Seat Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28" \$73.61
		\$85.87
		\$98.13

Common Cushion Upgrades			
Positioning:	<input type="radio"/> SR E0956	Lateral Support Reduction Cushions - Pair	\$114.64
	<input type="radio"/> SAS E0956	Side Arm Supports - Pair	\$151.37
	<input type="radio"/> CC-LTS E0956	Comfort Company Lateral Supports - Pair	\$540.92
	<input type="radio"/> LTS E0956	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED) - Pair	\$664.46

Front Mount Configuration	Wheel Options	Rear Mount Configuration	
<input type="radio"/> 5" Total Lock Caster Wheels - Pair	Standard	<input type="radio"/> 20" Mag Wheels - Pair Standard	
<input type="radio"/> 5" Locking Twin Caster Wheels - Pair	\$218.15	<input type="radio"/> 14" Mag Wheels - Pair N/C	
<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels - Pair*	\$218.15	<input type="radio"/> 22" Mag Wheels - Pair* N/C	
*(available with rear mags only)		*(22" mag may slightly limit tilt depending on arm support height)	
<input checked="" type="radio"/> Anti-Tip Devices(factory required safety feature) (E0971)	\$114.64	<input type="radio"/> Parking Brake Extension Handles (E0961)- Pair	\$73.46

Lower Leg Support Options (Select One)	
<input type="radio"/> Angle Adjustable, Articulating Split Lower Leg Support (E0990) (not available in 16" seat width)	\$247.09
Choose Number of Straps (foot support length range) <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")	\$247.09
<input type="radio"/> Angle Adjustable, Articulating Single Lower Leg Support (E0990)	\$247.09
Choose Number of Straps (foot support length range) <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")	
<input type="radio"/> Angle Adjustable, Articulating ABS Lower Leg Support (E0990) (includes foot cushions 16"-24" widths)	\$394.36

Foot Cushions: <input type="radio"/> FB-LONG E0994 Lower Leg Support Footbox - Full Extended Cushion <input type="radio"/> Omit Foot Cushions	\$421.98
<input type="radio"/> APP Calf & Foot APP Calf and Foot Support	Seat Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" \$188.91
	<input type="radio"/> 22" <input type="radio"/> 24" \$220.80
	<input type="radio"/> 26" <input type="radio"/> 28" \$252.69

Arm Support Options (Select One)	
Armrest Options: <input type="radio"/> Arm Supports - Pair	Standard
<input type="radio"/> Long Arm Supports - Pair	\$91.27
<input type="radio"/> Narrow Arm Support - Pair	\$91.27
<input type="radio"/> Cushioned Flip Arm Supports - Pair	\$140.00
<input type="radio"/> Long Cushioned Flip Arm Supports	\$180.00
<input type="radio"/> Arm Support Trough - Each	\$166.82
<input type="radio"/> Long Cushioned Arm Supports* - Pair	\$146.92
*Not compatible with a tray	
Arm Support Style: <input checked="" type="radio"/> Seat Mounted Adjustable Removable Arm Supports (E0973)	\$333.90

Common Accessories	
Pelvic Belts: <input type="radio"/> S-BELT E0978 Pelvic Belt <input type="radio"/> Center Buckle <input type="radio"/> Side Buckle	\$60.10
<input type="radio"/> S-BELT-A E0978 Pelvic Belt Auto Buckle <input type="radio"/> Center Buckle <input type="radio"/> Side Buckle	\$60.10
<input type="radio"/> S-BELT-PD E0978 Pelvic Belt Cushioned	\$86.81
<input type="radio"/> S-BELT-PD-A E0978 Pelvic Belt Cushioned & Auto Buckle	\$86.81
Thigh Belts: <input type="radio"/> T-BELT E0957 Thigh Belt Cushioned	\$140.24
<input type="radio"/> T-BELT-A E0957 Thigh Belt Cushioned & Auto Buckle	\$144.69
<input type="radio"/> T-BELT-NEO E0957 Thigh Belt Cushioned Neoprene	\$140.24
<input type="radio"/> T-BELT-NEO-AE0957 Thigh Belt Cushioned Neoprene Auto Buckle	\$144.69
<input type="radio"/> T-BELT-NEO-APP E0957 Thigh Belt Cushioned Neoprene for APP & SPP Packages	\$140.24
<input type="radio"/> T-BELT-NEO-A-APP E0957 Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages	\$144.94
Trays: <input type="radio"/> TRAY-ABS E0950 Tray ABS with Strap <input type="radio"/> Sewn Tray Strap <input type="radio"/> Grey <input type="radio"/> Black	\$357.27
<input type="radio"/> TRAY-ABS-OVR E0950 Tray ABS with Strap (24" - 28" widths) <input type="radio"/> Black Only	\$513.09
<input type="radio"/> TRAY-CLEAR E0950 Tray Clear with Strap <input type="radio"/> Sewn Tray Strap	\$369.52
<input type="radio"/> TRAY-HALF E0950 Half Tray (side specific)(includes upper side pad)	\$490.83
<input type="radio"/> TRAY-COVER Cushioned ABS Tray Cover	\$181.42
Stabilizers: <input type="radio"/> STABIL-FR <input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one)(pair)	\$103.51
<input type="radio"/> STABIL-LAT Lateral Frame Stabilizers	\$181.42
Other: <input type="radio"/> IV-POLE K0105 IV Pole and Assembly - Front Mount	\$114.64
<input type="radio"/> O2-HOLDER E2208 Oxygen Holder and Assembly - Front Mount	\$192.55
<input type="radio"/> ACCESSORY-BAG Accessory Bag (Black) Each	\$45.86
<input type="radio"/> BAG PERSONALIZING (2 Initials) <input type="radio"/> Brush For OR <input type="radio"/> Block Font	\$28.94
<input type="radio"/> O2-HOLDER-BAG (E2208) Oxygen Holder, Soft Carry Bag	\$97.52
<input type="radio"/> CUP-HOLDER Cup Holder	\$51.20

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorial services needed to deliver are requested when an order is submitted to Broda.

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.