

PO: _____ **Date:** _____ **Customer #:** _____

Customer Name: _____ Ship To Name: _____
 Address: _____ Address: _____
 Phone: _____ Phone: _____
 Contact: _____ Contact: _____
 Email: _____ Tag: _____



(shown with optional upgrades)

Standard Features	
• 250lb Maximum Occupant Weight (rocking)	• Up to 40° Infinitely Adjustable Seat Tilt (non-rocking)
• 350lb Maximum Occupant Weight (non-rocking)	• Rocking Mechanism (for up to 33° of tilt)
• 500 lb Maximum Occupant Weight with the Heavy Duty Package (non-rocking only)	• Up to 33° Infinitely Adjustable Seat Tilt (rocking)
• Tubular Steel, Powder Coated Frame	• Seat Mount Lower Leg Support Assembly
• Adjustable Seat Surface Height	• 2 x 20" Front Mag wheels with Front Anti-Tippers
• Independent Adjustable Arm Support Height	• 2 x 5" Total Lock Rear Caster Wheels

Encore Pedal Wheelchair® - Model: 48V4/48RV4 MSRP \$3,874.95

- WC19 Transit Tie Down Vehicle Package*** - Factory Installed Only \$375.13
*Available up to a 24" seat width, unless selected with the Heavy Duty Package.
- Heavy Duty Package**** - Factory Installed Only, Increases Maximum Occupant Weight to **500 lbs.** \$500.00
**Available in seat widths 22"-28" only. Non-rocking only. Aluminum front casters not available.
 Sub-mount seat hanger style only with pin & saddle, angle-adjustable footplates (F-SW-ANGLE) only.
 If combined with WC19 package, only mag wheels are available (20", 22" front or 14" rear with twin front casters)

Seat and Back Support Options

Rocking Mechanism:	<input type="radio"/> Eliminate Rocking Function -\$387.50	Seat Depth:	<input type="radio"/> 17" Standard
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C		<input type="radio"/> 18.5" N/C
	<input type="radio"/> 22" \$512.73		<input type="radio"/> 20" \$122.66
Non-Rocking Only	<input type="radio"/> 24" <input type="radio"/> 26" \$1,161.62	Standard Back	<input type="radio"/> 32" Standard
	<input type="radio"/> 28" \$1,161.62	Support Height:	<input type="radio"/> 33.5" \$122.66
Seat Surface	<input type="radio"/> 15" Standard		<input type="radio"/> 35" \$122.66
Height Rocking	<input type="radio"/> 14" <input type="radio"/> 16" <input type="radio"/> 17" N/C	Short Back	<input type="radio"/> 22" (does not include -\$490.67
	<input type="radio"/> 18" <input type="radio"/> 19" N/C	Support Height	Adjustable Lateral Shoulder Supports and Assembly, for Neck Support Option - see page 2)
Seat Surface	<input type="radio"/> 12.5" <input type="radio"/> 13.5" <input type="radio"/> 14.5" N/C		
Height No Rock	<input type="radio"/> 15.5" <input type="radio"/> 16.5" <input type="radio"/> 17.5" N/C		

Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair Standard	
<input type="radio"/> 22" Mag Wheels - Pair N/C	<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels* \$233.60	

*(available with front mags only.)

- 5" Total Lock Caster Wheels - Pair N/C
- 5" Total Lock Caster Wheels - Pair (\$680.01)
- 5" Fixed Non-Locking Caster Wheels - Pair N/C
- Parking Brake Extension Handles (pair) \$78.51
- Front-Mount Wheel Locks - Pair N/C

Arm Support Options

Arm Support	<input type="radio"/> Arm Support - Pair Standard	Arm Support	<input type="radio"/> Height Adjustable - Pair Standard
Cushion	<input type="radio"/> Arm Support Trough - Each \$166.82	Style:	<input type="radio"/> Height and Depth Adjustable - Pair \$169.27
Options:	<input type="radio"/> Cushioned Flip Arm Supports - Pair \$140.00	Arm Mount	<input checked="" type="radio"/> Seat Mounted Arm Supports - Pair Standard
	<input type="radio"/> Cushioned Arm Support* - Pair \$98.13		

**Not compatible with trays*

Seating and Frame Color Options			
Comfort Tension <input checked="" type="radio"/> Black	Standard	N/C	Frame Colors: <input type="radio"/> Black <input type="radio"/> Deep Orange <input type="radio"/> Tungsten Sand <input type="radio"/> Bronze <input type="radio"/> Zinc Yellow <input type="radio"/> Royal Purple <input type="radio"/> Carmine Red <input type="radio"/> Ocean Blue <input type="radio"/> Satin Nickel <input type="radio"/> Passion Pink
Seating Color:			
Upholstery Style Options (Select one)			
Accent Color Vinyl:	<input type="radio"/> Interlace Indigo***	<input type="radio"/> Interlace Espresso***	<input type="radio"/> Interlace Graphite***
	<input type="radio"/> Modena Black		
Vinyl:	<input type="radio"/> Beeline Portobello	<input type="radio"/> Beeline Geyser	<input type="radio"/> Beeline Forge <input type="radio"/> Kaleidoscope Smokehouse
Nylon:	<input type="radio"/> Retro Camo Jungle		\$159.47
Upholstery Colors come with seat and back cushions in black vinyl. *Interlace Indigo, Espresso, and Graphite available with seat and back cushions in matching vinyl color			
Cushions - Choose Option 1, 2, OR 3			
<input type="radio"/> Option 1 - Broda Standard Cushions			N/C
Includes: Neck Support, Adjustable Lateral Shoulder Supports, Lower Side Covers, Seat and Back Support Cushions with Covers			
Seat & Back Pad Cover	<input type="radio"/> Black Vinyl		Standard
Options:	<input type="radio"/> Accent Color Vinyl ***Option For Interlace: Indigo, Espresso or Graphite only		N/C
	(Select Black Vinyl OR Accent Color Vinyl Above)		
	<input type="radio"/> Full Vinyl Wrap Seat and Back Support Cushions		\$122.66
	(Select Black Vinyl OR Accent Color Vinyl Above)		
	<input type="radio"/> Grey Terry Cover		N/C
Lateral Shoulder Support	<input type="radio"/> Lateral Shoulder Supports - Pair		N/C
Options:	<input type="radio"/> Long Lateral Shoulder Supports - Pair		\$166.81
Optional Substitutions Available Between APP & SPP Package Cushion Components			
<input type="radio"/> Option 2 - Broda Additional Positioning Padding (APP), 24"-28" seat widths for non-rocking only			
Broda APP Package Includes: Seat, Back and Head Support Cushions (no Head Support or Lateral Shoulder Supports on short back), plus Adjustable, Removable Lateral Shoulder Supports and Lateral Upper Leg Support Cushions			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21" <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21" <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"
	<input type="radio"/> 22" <input type="radio"/> 24"		\$451.41
	<input type="radio"/> 26" <input type="radio"/> 28"		\$706.54
			\$836.57
<input type="radio"/> Option 3 - Broda Slim Posture Padding (SPP), 24"-26" seat widths for non-rocking only			
Broda SPP Cushion Options: Seat, Back Support, and Back Support Extension Cushions			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20" <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20" <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"
	<input type="radio"/> 22" <input type="radio"/> 24"		\$161.97
	<input type="radio"/> 26"		\$178.64
Lower Back Support	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Upper Back Support	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24"
	\$142.91		\$107.19
	\$171.49		\$119.09
13" Height	<input type="radio"/> 26"	12" Height	<input type="radio"/> 26"
	\$207.22		\$131.01
Matrx® Head Support (Available with 22" Short Back Option Only)			
Head Support Mounting:	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support		\$235.40
Choose One	<input type="radio"/> Elan Flipp Head Support Hardware - Required with Selection of Head Support		\$535.00
	<input type="radio"/> Loxx Head Support Hardware (E0128) - Required with Selection of Head Support		\$215.00
	<input type="radio"/> Loxx Flipp Head Support Hardware (E0128) - Required with Selection of Head Support		\$515
Support Options:	<input type="radio"/> Head Support Cushion	<input type="radio"/> 10" <input type="radio"/> 14"	\$235.00
	<input type="radio"/> Hug Head Support Cushion	<input type="radio"/> 10" <input type="radio"/> 14"	\$235.00
	<input type="radio"/> Head Support 4 Point Cushion	<input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10"	\$325.00



Encore Pedal Wheelchair Model: 48V4/48RV4

US MSRP Order Form
April 2024

Headband with 4 Point Cushion

Small

Large

\$450.00

Drop In Style	Must Choose One Lower Leg Support Assembly Style	Pin and Saddle Style
<input type="radio"/> Lower Leg Support Assembly (Drop-in)- Pair N/C	<input type="radio"/> Lower Leg Support Assembly (Pin & Saddle) - Pair* N/C	

*Pin & Saddle Lower Leg Support Assembly not available with 16" SW chair due to overlapping foot supports.

Must Choose One Seat Mount Lower Leg Support Assembly Location	
<input type="radio"/> Seat Frame Mount* N/C	<input type="radio"/> Sub-Seat Frame Mount** N/C

*Not available on rocking version with 17" seat depth

**Improves the ease of the rocking mechanism


*Not available on non-rocking version when choosing 12.5"

seat surface height with pin & saddle style lower leg support assembly

Seat Frame Mount



Sub-Seat Frame Mount



Drop-In Style	Lower Leg Support Options	Pin & Saddle Style
<input type="radio"/> F-SW-IKN: Swing Away Lower Leg Supports Silver - Pair \$116.54	<input type="radio"/> F-SW: Swing Away Lower Leg Supports Black- Pair \$186.45	
<input type="radio"/> F-SW-ELV-IKN: Elevating Swing Away Lower Leg Supports \$484.53 Silver Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair	<input type="radio"/> F-SW-ELV: Elevating Swing Away Lower Leg Supports \$250.23 Black w/ Flip-up Foot Supports - Pair	
<input type="radio"/> F-SW-IKN-ADJ: Swing Away Lower Leg Supports Black \$239.37 Steel w/ Adjustable Flip-up Foot Supports - Pair	<input type="radio"/> F-SW-ANGLE Swing Away 70 Degree Lower Le \$316.48 Support w/ Angle Adjustable Footplates Black*	
	<input type="radio"/> LEG-FOOT-ENCORE*: Lower Leg Support Foot Board \$291.95	
	*Available with F-SW & F-SW-ELV only	
	<input type="radio"/> LEGPANEL-ENCORE: Lower Leg Support ABS Cushions \$98.13	
<input type="radio"/> 2" Extended Footrest Tubes - <i>Only compatible on the F-SW-IKN and F-SW-IKN-ADJ</i>		\$71.45

Encore Pedal with Seat Mounted Hangers	Compatible with 14.5" Seat Height and above
Encore Pedal with Sub-Seat Mounted Hangers	Compatible with 15.5" Seat Height and above
Encore Pedal with Rocker with Seat Mounted Hangers	Compatible with 15" Seat Height and above
Encore Pedal with Rocker with Sub-Seat Mounted Hangers	Compatible with 16" Seat Height and above

Common Accessories				
Head	<input type="radio"/> NR-CV	Head Support Visco Foam Contoured Cushion	<input type="radio"/> Triptex	<input type="radio"/> Omit Standard Head Support \$130.03
Supports:	<input type="radio"/> NS	Neck Support Soft Visco Cushion		\$130.03
Positioning	<input type="radio"/> SR	Width Reduction Cushion- Pair		\$122.66
Supports:	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports - Pair		\$578.98
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)		\$711.45
Pelvic Belts:	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle \$63.79
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> CenterBuckle	<input type="radio"/> Side Buckle \$63.79
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned		\$93.23
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned and Auto Buckle		\$93.23
Thigh Belts:	<input type="radio"/> T-BELT	Thigh Belt Cushioned		\$149.65
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned and Auto Buckle		\$154.82
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene		\$149.65
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene and Auto Buckle		\$154.82
	<input type="radio"/> T-BELT-NEO-APP	Thigh Belt Cushioned Neoprene for APP & SPP Packages		\$149.65
	<input type="radio"/> T-BELT-NEO-A-API	Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages		\$154.82

Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$382.71
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)		<input type="radio"/> Black Only	\$549.55
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$394.98
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes lower side cover)		<input type="radio"/> Grey <input type="radio"/> Black	\$525.01
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover			\$193.82
Other:	<input type="radio"/> IV-POLE*	IV Pole & Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount		\$122.66
	<input type="radio"/> O2-HOLDER*	Oxygen Holder & Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount		\$206.03
	*If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted.				
	<input type="radio"/> STABIL-FR	Rear Stabilizer Frame			\$110.40
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each			\$49.70
	<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font	<input type="radio"/> Block Font		\$30.68
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag		(Not recommended with 22" back height due to potential ground interference.)	\$97.52
	<input type="radio"/> OCUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)			\$55.20

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorial services needed to deliver are requested when an order is submitted to Broda.

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.