

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
• Comfort Tension Seating®	• 350lb Maximum Occupant Weight
• 2 x 20" Rear Mag Wheels	• 2 x 5" Total Lock Caster Wheels
• Up to 38° Infinitely Adjustable Seat Tilt	• Adjustable Seat Surface Height
• Up to 90° Infinitely Adjustable Back Support Recline	• Height Adjustable Arm Supports
• Tubular Steel, Powder Coated Frame	• Angle Adjustable, Articulating Lower Leg Support

Midline Positioning Wheelchair - Model: MID MSRP \$4,300.59

Seat and Back Support Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C	Seat Depth:	<input type="radio"/> 15.5" \$122.66
	<input type="radio"/> 22" \$512.73		<input type="radio"/> 17" Standard
	<input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28" \$1,160.40		<input type="radio"/> 18.5" N/C
Seat Surface Height:	<input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" N/C		<input type="radio"/> 20" \$122.66
	<input type="radio"/> 20" Standard	Back Support Height:	<input type="radio"/> 32" Standard
	<input type="radio"/> 21" <input type="radio"/> 22" N/C		<input type="radio"/> 33.5" \$122.66
	*Seat heights below 20" (51 cm) reduce the available range of tilt.		<input type="radio"/> 35" \$122.66

Front Options	Wheel Options	Rear Options
<input type="radio"/> 5" Total Lock Caster Wheels - Pair Standard	<input type="radio"/> 20" Mag Wheels - Pair Standard	
<input type="radio"/> 20" Mag Wheels - Pair N/C	<input type="radio"/> 1x 5" Directional Locking 1x Total Lock Caster Wheel (\$680.01)	
<input type="radio"/> 22" Mag Wheels - Pair N/C	<input type="radio"/> 22" Mag Wheels - Pair N/C	
<input type="radio"/> 24" Mag Wheels - Pair N/C	<input type="radio"/> Parking Brake Extension Handles - Pair \$78.51	

Lower Leg Support Options (Select One)			
<input type="radio"/> Angle Adjustable, Articulating Split Lower Leg Support (not available in 16" seat width)			Standard
Choose Number of Straps (support length range)	<input type="radio"/> 2 (12"-19")	<input type="radio"/> 3 (14"-21")	<input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Angle Adjustable, Articulating Single Lower Leg Support			N/C
Choose Number of Straps (support length range)	<input type="radio"/> 2 (12"-19")	<input type="radio"/> 3 (14"-21")	<input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")

Arm Support Options			
Arm Support	<input type="radio"/> Arm Supports - Pair Standard	Arm Support	<input type="radio"/> Height Adjustable Arm Supports - Pair Standard
Cushion	<input type="radio"/> Arm Support Trough - Each \$166.82	Styles:	<input type="radio"/> Height & Depth Adjustable Arm Supports Pair \$169.
Options:	<input type="radio"/> Cushioned Flip Arm Supports - Pair \$140.00	Arm Support	<input checked="" type="radio"/> Base Mounted Arm Supports Standard
	<input type="radio"/> Cushioned Arm Supports* -Pair \$98.13	Mount:	
	*Not compatible with a tray		

Seating and Frame Color Options			
Comfort Tension	<input type="radio"/> Grey <input type="radio"/> Blue N/C	Cushion	<input type="radio"/> Burgundy <input type="radio"/> Charcoal N/C
Seating Color:	<input type="radio"/> Black <input type="radio"/> Green N/C	Color:	<input type="radio"/> Sapphire Blue N/C
Frame Color:	<input type="radio"/> Grey <input type="radio"/> Black N/C		

Cushions- Choose Option 1, 2, OR 3

Option 1 - Broda Cushions N/C

Standard Broda Cushions:

Includes: Neck Support, Lateral Shoulder Supports, Foot Support Cushions, Upper & Lower Side Covers, Seat & Back Cushions with Covers

Seat & Back Cushion Grey Terry Cover Startex Cover Vinyl Cover N/C

Cover Options: Full Vinyl Wrap Seat & Back Cushions \$122.66

Lateral Shoulder Support Standard Short Rounded Short ACC Short N/C

Options: Standard Long Rounded Long ACC Long \$166.82

Huntington's Special Includes Long Lateral Shoulder Support Covers, Thigh Belt & Long Lower Leg Support Footbc \$667.30

Padding (HSP) package:

Optional Substitutions Available Between APP & SPP Components

Option 2 - Broda Additional Positioning Padding (APP) (not compatible with 35" back support height)

Broda APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Covers

Seat Width: 16" 18" 20" **Seat Depth:** 17" 19" 21" \$640.31

22" 24" 17" 19" 21" \$770.33

26" 28" 17" 19" 21" \$900.36

Option 3 - Broda Slim Posture Padding (SPP)

Broda SPP Padding Options: Seat, Back Support, and Back Support Extension Cushions

Seat Width: 16" 18" 20" **Seat Depth:** 17" 18.5" 20" \$161.97

22" 24" 17" 18.5" 20" \$178.64

26" 17" 18.5" 20" \$207.22

Back Support 16" 18" 20" \$142.91 **Back Support** 16" 18" 20" \$107.19

13" Height 22" 24" \$171.49 **Extension:** 22" 24" \$119.09

26" \$207.22 12" Height 26" \$131.01

Acute Care Configuration (ACC):

Includes: Full Vinyl Covers, Acute Care Shoulder Supports, ABS Lower Leg Support with Foot Support Cushions \$583.88

(Factory installed)(skip lower leg support cushion package section above when selecting this option)(16"-24" widths)

Common Accessories

Neck Supports: NR-CV Neck Support Visco Foam Contoured Cushion Triptex Omit Standard \$130.03

NS Neck Support Soft Visco Cushion Neck Support \$130.03

Positioning: SR Width Reduction Cushions - Pair \$122.66

SAS Side Arm Supports - Pair \$161.92

CC-LTS Comfort Company Lateral Supports - Pair \$578.98

LTS Lateral Torso Supports *(FACTORY INSTALLED RECOMMENDED)* - Pair \$711.45

Foot Support: FB-LONG Lower Leg Support Footbox - Full Extended Cushion Omit Foot Support \$421.98

Pelvic Belts: S-BELT Pelvic Belt Center Buckle Side Buckle \$63.79

S-BELT-A Pelvic Belt Auto Buckle Center Buckle Side Buckle \$63.79

S-BELT-PD Pelvic Belt Cushioned \$93.23

S-BELT-PD-A Pelvic Belt Cushioned & Auto Buckle \$93.23

Thigh Belts: T-BELT Thigh Belt Cushioned \$149.65

T-BELT-A Thigh Belt Cushioned & Auto Buckle \$154.82

T-BELT-NEO Thigh Belt Cushioned Neoprene \$149.65

T-BELT-NEO-A Thigh Belt Cushioned Neoprene Auto Buckle \$154.82

T-BELT-NEO-APP Thigh Belt Cushioned Neoprene for APP & SPP Packages \$149.65

T-BELT-NEO-A-API Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages \$154.82

Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey	<input type="radio"/> Black	\$382.71
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)		<input type="radio"/> Black Only		\$549.55
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap			\$394.98
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes upper side pad)				\$525.01
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover				\$193.82
Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one)(pair)				\$110.40
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers				\$193.82
Other:	<input type="radio"/> IV-POLE*	IV Pole and Assembly	<input type="radio"/> Front	or	<input type="radio"/> Rear Mount	\$122.66
	<input type="radio"/> O2-HOLDER*	Oxygen Holder and Assembly	<input type="radio"/> Front	or	<input type="radio"/> Rear Mount	\$206.03
	*IV-POLE and O2-HOLDER not compatible with rear mounted mag wheels					
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each				\$49.70
	<input type="radio"/> BAG PERSONALIZING	(2 Initials: <input type="radio"/> Brush Fon OR <input type="radio"/> Block Font)				\$30.68
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag				\$97.52
	<input type="radio"/> CUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)				\$55.20

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
 - Residential Delivery Service:** Delivery to a residential location
 - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
 - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
 - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**