

**PO:** \_\_\_\_\_ **Date:** \_\_\_\_\_ **Customer #:** \_\_\_\_\_

Customer Name: \_\_\_\_\_ Ship To Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Contact: \_\_\_\_\_  
 Email: \_\_\_\_\_ Tag: \_\_\_\_\_



Standard Features	
• Up to 14° Infinitely Adjustable Seat Tilt (Posterior)	• Up to 800lb Maximum Occupant Weight
• Up to 8° Infinitely Adjustable Seat Tilt (Anterior)	• Contoured Backrest
• Up to 65° Infinitely Adjustable Back Support Recline	• Comfort Tension Seating®
• Stainless Steel, Powder Coated Frame	• Width Adjusting Removable Arm Supports
• 4 x 6" Total Lock Stainless Steel Front & Rear Caster Wheels	• Head Support
• 2 x 8" Locking Middle Wheels with One-Touch Emergency Brake	

**Revive Bariatric Shower Commode - Model: B385 MSRP \$7,958.43**  
**Preset Factory Installed Cylinder Lift Forces(Select one)**

Cylinder Configuration	Suggested Occupant Weight (lbs)	Cylinder Lift Force	
2 x 300N	<input type="radio"/> 100-350 lbs.	135 lbs.	N/C
2 x 600N	<input type="radio"/> 350-500 lbs.	270 lbs.	N/C
2 x 800N	<input type="radio"/> 500-650 lbs.	360 lbs.	N/C
2 x 300N + 2 x 600N	<input type="radio"/> 650-800 lbs.	405 lbs.	N/C

\*Guideline Only. Contact Broda for additional weight accomodation options

Seat and Back Support Options			
<b>Seat Width:</b> <input checked="" type="radio"/> 27" Standard	<b>Seat Depth:</b> <input type="radio"/> 21" N/C		
(33" with width adjustable armrests)	<input type="radio"/> 23" N/C		
<b>Seat Surface:</b> <input type="radio"/> 21" N/C	<b>Back Height:</b> <input checked="" type="radio"/> 26" Standard		
<b>Height:</b> <input type="radio"/> 23" N/C	<b>No Anterior Tilt*:</b> <input type="radio"/> N/C		

\*Max Upright Position - 5 degree posterior tilt

Wheel Options	
<input checked="" type="radio"/> 6" Total Lock Stainless Steel Front & Rear Caster Wheels & 8" Locking Middle Wheels with One-Touch Emergency Brake	Standard

Seat Options (Select one)	Bedpan & Basket Options
<input type="radio"/> Seat Commode Cushioned Standard	<input checked="" type="radio"/> Bedpan with Basket Standard
<input type="radio"/> Seat Transport Cushioned N/C	

Arm Support Options	
<b>Arm Support Options:</b> <input checked="" type="radio"/> Arm Supports	Standard
<b>Arm Support Style:</b> <input checked="" type="radio"/> Seat Width Extending (extends width to 33")	Standard

Head Support	
<b>Head Support:</b> <input checked="" type="radio"/> Molded 2" Offset	Standard

Strap Color	
<b>Strap Color:</b> <input type="radio"/> Grey <input type="radio"/> Blue <input type="radio"/> Green	N/C

Common Accessories	
<b>Lateral Supports:</b> <input type="radio"/> Molded Laterals (pair)	\$98.13
<b>Pelvic Belt:</b> <input type="radio"/> Pelvic Belt (side buckle)(large size)	\$63.79

**Accessorial Services:** Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

**any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

**Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**

---