

**PO:** \_\_\_\_\_ **Date:** \_\_\_\_\_ **Customer #:** \_\_\_\_\_

Customer Name: \_\_\_\_\_ Ship To Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Contact: \_\_\_\_\_  
 Email: \_\_\_\_\_ Tag: \_\_\_\_\_



Standard Features	
• Height Adjustable, Seat or Base Mount Arm Supports	• 600lb Maximum Occupant Weight
• Tubular Steel, Powder Coated Frame	• Foot Pump Operated Tilt
• 2 x 8" Rear Directional Locking Caster Wheels	• 2 x 6" Total Lock Front Caster Wheels
• Up to 30° Infinitely Adjustable Seat Tilt	• Stand Assist Handles
• Up to 45° Infinitely Adjustable Back Support Recline	• Comfort Tension Seating®
• Up to 10° Infinitely Adjustable Anterior Seat Tilt	

**Vanguard® Bariatric Wheelchair - Model: 985 VG MSRP \$7,764.62**

**Seat and Back Support Options**

<b>Seat Width:</b>	<input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28" <input type="radio"/> 30"	N/C	<b>Seat Depth:</b>	<input type="radio"/> 21" <input type="radio"/> 23"	Standard N/C
<b>Seat Surface</b>	<input type="radio"/> 20"	N/C	<b>Back Support</b>	<input type="radio"/> 32" <input type="radio"/> 34"	N/C
<b>Height:</b>	<input type="radio"/> 21" <input type="radio"/> 22" <input type="radio"/> 23" <input type="radio"/> 24"	Standard N/C N/C N/C	<b>Height:</b>	<input type="radio"/> 36" <input type="radio"/> 38" <input type="radio"/> 40"	N/C N/C N/C

Back Option		
<input type="radio"/> Straight Back	N/C <input type="radio"/> Curved Back ( <i>concave</i> )	N/C

Wheel Options	
<input checked="" type="radio"/> 8" Multi-Directional Locking Rear Caster Wheels & 6" Total Lock Front Caster Wheels	Standard

Lower Leg Support Options ( <i>Select one</i> )	
<input checked="" type="radio"/> Angle Adjustable, Articulating Single Lower Leg Support	Standard

Arm Support Options		
<b>Arm Support</b>	<input type="radio"/> Arm Support - Pair Standard	<b>Arm Support:</b> <input checked="" type="radio"/> Height Adjustable Arm Support Standard
<b>Options:</b>	<input type="radio"/> Arm Support Trough - Each \$166.82	<b>Arm Support</b> <input type="radio"/> Seat Mounted Arm Support Standard
	<input type="radio"/> Cushioned Flip Arm Supports \$140.00	<b>Mount:</b> <input type="radio"/> Base Mounted Arm Support Standard
	<input type="radio"/> Long Cushioned Flip Arm Supports \$180.00	
	<input type="radio"/> Cushioned Arm Supports* -Pair \$98.13	

*\*Not compatible with a tray*

Seating and Frame Color Options			
<b>Strap Color:</b>	<input type="radio"/> Grey <input type="radio"/> Blue N/C	<b>Cushion</b>	<input type="radio"/> Burgundy <input type="radio"/> Charcoal N/C
<b>Frame Color:</b>	<input checked="" type="radio"/> Grey Standard	<b>Color:</b>	<input type="radio"/> Sapphire Blue N/C

Cushions	
<input type="radio"/> Standard Broda Cushions:	N/C
Includes: Neck Support, Lateral Shoulder Supports, Calf Support Cushion, Upper & Lower Side Covers, Seat & Back Cushions with Covers	
<input type="radio"/> Acute Care Configuration (ACC): Full Vinyl Wrap Seat and Back Cushions	\$305.44
<b>Seat, Back Cushion Cover</b>	<input type="radio"/> Grey Terry Cover <input type="radio"/> Startex Cover N/C
<b>Options:</b>	<input type="radio"/> Vinyl N/C

Common Accessories				
<b>Pelvic Belts:</b>	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle <input type="radio"/> Side Buckle	\$63.79
(26" - 28" width only)	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> Center Buckle <input type="radio"/> Side Buckle	\$63.79
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned		\$93.23
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned & Auto Buckle		\$93.23
<b>Thigh Belts:</b>	<input type="radio"/> T-BELT	Thigh Belt Cushioned		\$149.65
(26" - 28" width only)	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned & Auto Buckle		\$154.82
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene		\$149.65
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene Auto Buckle		\$154.82
<b>Trays:</b>	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 30" width)	<input type="radio"/> Black Only	\$549.55
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes upper side pad)		\$525.01
	<input type="radio"/> TRAY-COVER	Cushioned ABS Tray Cover (26" or 28" width only)		\$193.82
<b>Other:</b>	<input type="radio"/> IV-POLE	IV Pole and Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount	\$122.66
	<input type="radio"/> O2-HOLDER	Oxygen Holder and Assembly	<input type="radio"/> Front or <input type="radio"/> Rear Mount	\$206.03
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each		\$49.70
	<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font OR <input type="radio"/> Block Font		\$30.68
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag		\$97.52
	<input type="radio"/> CUP-HOLDER	Cup Holder		\$55.20

Common Cushion Upgrades				
<b>Neck Support:</b>	<input type="radio"/> NS	Neck Support (one size)	<input type="radio"/> Omit Standard Neck Support	\$130.02
<b>Positioning:</b>	<input type="radio"/> SR	Width Reduction Cushions - Pair		\$114.64
<b>Foot Support:</b>	<input type="radio"/> SP-SINGLE-985	Sole Single 985 Cushion		\$98.13

**Additional Options:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

**Delivery Instructions or Special Requests:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Total: \$  

**Accessorial Services:** Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

**any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

**Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**