

## Patient, Facility, & Caregiver Benefits

- Tilt and recline provides ease of use for frequent repositioning
- Design provides safety and stability helping to reduce risk of injury
- Proper seating alignment decreases postural deviations reducing the use of restraints




Posterior tilt aids in pressure redistribution



Tilt with semi-recline accommodates hip flexion limitations, allowing patients to be out of bed



Optional contracture sling accommodates various levels of knee contractures



Comfort Tension Seating® improves postural support and sitting tolerance

### Essential Positioning & Peace of Mind

#### SPECIFICATIONS

<b>INFINITELY ADJUSTABLE TILT</b>	Up to 25°						
<b>DEGREES OF RECLINE</b>	Up to 45°						
<b>WEIGHT CAPACITY</b>	250 lbs. (114 kg)						
<b>SEAT WIDTHS</b>	16 in (40 cm)	18 in (45 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)	28 in (71 cm)
<b>SEAT DEPTHS</b>	15.5 in (39 cm)	17 in (43 cm)	18.5 in (45 cm)	20 in (58 cm)			
<b>SEAT SURFACE HEIGHTS (ADJUSTABLE)</b>	19.5 in (49 cm)	20.5 in (51 cm)	21.5 in (54 cm)				
<b>BACK HEIGHT</b>	32 in (82 cm)	33.5 in (86 cm)	35 in (89 cm)				
<b>WHEELS - CASTER WHEELS</b>	<b>FRONT or REAR</b> Total Lock Caster Wheels 5 in (12.5 cm)						
	<b>DIRECTIONAL</b> Rear Caster Wheels 8 in (20 cm)						
<b>WHEELS - MAG WHEELS</b>	<b>FRONT or REAR</b> Mount Mag Wheels 20 in (50 cm), 22 in (55 cm), 24 in (60cm)						
<b>LIFETIME WARRANTY FRAME</b>	16-Gauge Tubular Steel Frame <b>10 YEARS</b>						
<b>PRESSURE-REDUCING BACK SUPPORT</b>	Integrated with Comfort Tension Seating®						



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)

