

Patient, Facility, & Caregiver Benefits

- Pressure redistribution that aids in maintaining skin integrity
- Durable frame to withstand heavy workload healthcare environments
- Adaptable for wide range of patients and care needs



Infinitely-adjustable seat tilt maintains user comfort with smooth, effortless tilt.



Posterior tilt optimizes pressure redistribution.



Solid seat pan offers support and is compatible with a variety of seat cushion options.



Range of simple adjustments accommodates a variety of needs.

Shown with optional head support, **matrx**® Elite back support, **matrx**® Vi seat cushion, and angle-adjustable lower leg supports

Simplified Durability for Rehabilitation Care

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 45°				
ANGLE ADJUSTABLE BACK CANES	5°	10°	15°	20°	
WEIGHT CAPACITY	250 lbs. (114 kg)				
SEAT WIDTHS	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	
SEAT DEPTHS (ADJUSTABLE)	16 in (40 cm)	17 in (43 cm)	18 in (45 cm)	19 in (48 cm)	20 in (50 cm)
SEAT SURFACE HEIGHTS (ADJUSTABLE)	14 in (36 cm)	15 in (38 cm)	15.5 in (39 cm)	16 in (41 cm)	16.5 in (42 cm)
	17 in (43 cm)	17.5 in (44 cm)	18 in (45 cm)	19 in (48 cm)	
WHEELS - MAG WHEELS	REAR MOUNT 20 in (50 cm), 22 in (55 cm), or 24 in (60 cm)				
WHEELS - CASTER WHEELS	FRONT Caster Wheels with Fork 5 in (13 cm), 6 in (15 cm), or 8 in (20 cm)				
LIFETIME WARRANTY FRAME	16-Gauge tubular steel frame 10 YEARS				
GOVERNMENT CODING (US & CA)	HCPCS: E1161 ADP (ON): WMA5B-025				



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)



Matrix® seating options available to ship preconfigured on-chair.



LIBRA CUSHION

Foam
HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.



E2 BACK

Shell
Contoured aluminum with cutouts for lifting and breathability. Available in standard (3") OR deep (6") contours.



LOXX HEADREST

Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.