

Patient, Facility, & Caregiver Benefits

- Facilitates positioning for comfort and respiratory function
- Infinitely adjustable seat tilt provides superior pressure redistribution
- Proper seating alignment decreases postural deviations
- Optional recline feature adds even more positioning options for optimal user comfort




Tilt and recline functionality allow easy repositioning of users



Removable arm supports make transfers safer for caregivers and patients



Posterior Tilt effectively opens the diaphragm improving digestion, oxygenation, and organ function



Comfort Tension Seating[®] increases support and sitting tolerance

Enhanced Positioning & Functionality

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 40° Posterior						
INFINITELY ADJUSTABLE RECLINE	Up to 90°						
WEIGHT CAPACITY	350 lbs. (158 kg)						
SEAT WIDTHS	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)	28 in (65 cm)
SEAT DEPTHS	15.5 in (37.5 cm)	17 in (43 cm)	18.5 in (46 cm)	20 in (50 cm)	22 in (56 cm)		
SEAT SURFACE HEIGHTS	19 in (48 cm)	20.5 in (51 cm)	22 in (55 cm)	23.5 in (59 cm)			
WHEELS - MAG WHEELS	REAR MOUNT 14 in (36 cm), 20 in (51 cm)						
	MID MOUNT 14 in (36 cm)						
WHEELS - CASTER WHEELS	FRONT total lock caster wheels 5 in (13 cm)						
	REAR multi-directional locking caster wheels 6 in (15 cm)						
LIFETIME WARRANTY FRAME	16-Gauge tubular steel frame 10 YEARS						




Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)

Optional Accessories



O2 Holder



Lateral Support



Half Tray



Clear Tray



Padded Thigh Strap