

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive, stable positioning



Adjustable seat surface height goes as low as 12.5 in. (32 cm) for optimal heel strike



Front-pivot seat tilt maintains a proper position for self-propulsion



Optional Dynamic Rocking feature helps the user remain calm and content



Comfort Tension Seating® increases support as well as sitting tolerance

NEW

Heavy Duty Package
Increases weight limit to 500 lbs.

Mobility with Comfort, Functionality, & Style

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 40°							
WEIGHT CAPACITY	Up to 350 lbs. (159 kg)							
	Dynamic Rocking Up to 250 lbs. (114 kg)							
	Heavy Duty Package Up to 500 lbs. (227 kg)							
SEAT WIDTHS	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)	28 in (71 cm)	
	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	with Dynamic Rocking			
SEAT DEPTHS	17 in (43 cm)	18.5 in (46 cm)	20 in (50 cm)					
	15.5 in (39 cm)	17 in (50 cm)	18.5 in (46 cm)	20 in (50 cm)	with Dynamic Rocking			
BACK LENGTHS	22 in (56 cm)	32 in (81 cm)	33.5 in (85 cm)	35 in (89 cm)				
SEAT SURFACE HEIGHTS	ADJUSTABLE 12.5 in (31 cm) - 17.5 in (44 cm)							
	ADJUSTABLE 14 in (35 cm) - 19 in (48 cm) with Dynamic Rocking							
WHEELS - MAG WHEELS	FRONT 20 in (50 cm) or 22 in (55 cm)							
	REAR 14 in (36 cm)							
WHEELS - CASTER WHEELS	FRONT Total Lock Caster Wheels 5 in (12.5 cm)							
	REAR Casters 5 in (12.5 cm)							
LIFETIME FRAME WARRANTY	16-Gauge Tubular Steel Frame 10 YEARS							
PRESSURE-REDUCING SEAT & BACKREST	Integrated with Comfort Tension Seating®							



Low back option shown here with the **matrix** Elan head support.

Additional accessories available, including foot supports, trays, positioning straps, and more!

[Click here to see accessories, color options, order forms, and more!](#)