



## Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive, stable positioning



Shown with Libra cushion, E2 back, and Loxx Head support



**Back canes fold down** for easy storage options



**Dynamic rocking function** soothes and stimulates the user



Front-pivot seat tilt wmaintains a proper position for self-propulsion



**Open seating system** allows for almost any seating configuration





# Mobility with Comfort, Functionality, & Style

## **SPECIFICATIONS**

INFINITELY ADJUSTABLE TILT	Up to 33° (with Dynamic Rocking)				
WEIGHT CAPACITY	250 lbs. (114 kg)				
SEAT WIDTHS	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	
SEAT DEPTHS	16 in (40 cm)	17 in (43 cm)	18 in (45 cm)		
	19 in (48 cm)	20 in (50 cm)			
BACK LENGTHS	16 in (40 cm)	18 in (45 cm)	21 in (53 cm)	24 in (60 cm)	Angle adjustable to 5°, 10°, 15°, 20°
SEAT SURFACE HEIGHTS	ADJUSTABLE 14 in (35 cm) - 19 in (48 cm)				
WHEELS - MAG WHEELS	FRONT Mag Wheels 20 in (51 cm) or 22 in (55 cm)				
WHEELS - CASTER WHEELS	REAR Total Lock Caster Wheels 5 in (12.5 cm)				
	REAR Aluminum Non-Locking Wheels 5 in (12.5 cm)				
LIFETIME FRAME WARRANTY	16-Gauge Tubular Steel Frame 10 YEARS				
PRESSURE-REDUCING SEAT & BACKREST	Integrated with Comfort Tension Seating®				



Additional accessories available, including foot supports, trays, positioning straps, and more!

Click here to see accessories, color options, order forms, and more!



## LIBRA CUSHION

#### Foam

HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.



## **MULTI CUSHION**

**Three-part foam combination** Cooling gel foam and immersive support.



## LOXX HEADREST

Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.
Flipp hardware available.



## **E2 BACK**

**Shell** - Contoured aluminum with cutouts for lifting and breathability. Available in standard (3", 8 cm) OR deep (6", 15 cm) contours.

**Foam** - HR (high resilient) foam enclosed in cover.