

Patient, Facility, & Caregiver Benefits

- Caregiver activated tilt allows for easy, frequent repositioning
- Design provides safety and stability helping to reduce risk of injury
- Dynamic Rocking provides stimulation; excellent for memory care



Shown with optional rear mag wheels



Rocking motion helps calm agitation, creating a soothing effect



Posterior tilt aids in pressure redistribution



Front-pivot seat tilt enhances visual orientation and encourages self-mobility with proper foot-on-floor positioning



Comfort Tension Seating® improves postural support and sitting tolerance

Enhanced Mobility & Freedom for Memory Care

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 33°					
WEIGHT CAPACITY	250 lbs. (113 kg)					
SEAT WIDTHS	15.5 in (39 cm)	16 in (40 cm)	17 in (43 cm)	18.5 in (47 cm)	20 in (50 cm)	
SEAT DEPTHS	17 in (43 cm)	18.5 in (46 cm)	20 in (50 cm)			
SEAT SURFACE HEIGHTS	14 in (35 cm)	15 in (38 cm)	16 in (40 cm)	17 in (43 cm)	18 in (45 cm)	19 in (48 cm)
WHEELS - MAG WHEELS	REAR 20 in (35 cm) or 22 in (55 cm)					
WHEELS - CASTER WHEELS	FRONT Total Lock Caster Wheels 5 in (12 cm)					
	FIXED REAR Caster Wheels 5 in (12 cm)					
LIFETIME WARRANTY FRAME	16-Gauge Tubular Steel Frame 10 YEARS					
PRESSURE-REDUCING SEAT & BACK SUPPORT	Integrated with Comfort Tension Seating*					



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)

