BRODA®



Patient, Facility, & Caregiver Benefits

- Infinitely adjustable tilt makes it quick and simple to reposition the user
- Designed for safety, stability, and reduced risk of injury
- Provides stimulation and contentment excellent for memory care





Dynamic rocking feature helps calm agitation and absorb rigid movements



Posterior tilt improves posture and redistributes pressure



Front-pivot seat tilt enhances visual orientation, and encourages self-mobility with proper heel strike



Comfort Tension Seating[®] improves postural support and sitting tolerance

BRODA



Enhanced Mobility & Freedom for Memory Care

SPECIFICATIONS

| INFINITELY ADJUSTABLE TILT | Up to 33° (with dynamic rocking) | | | | | |
|-------------------------------|--|------------------|--------------------|------------------|------------------|------------------|
| WEIGHT CAPACITY | 250 lbs. (113 kg) | | | | | |
| SEAT WIDTHS | 16 in (40 cm) | 18 in (40 cm) | 20 in (50 cm) | 22 in (55 cm) | | |
| SEAT DEPTHS | 15.5 in (39cm) | 17 in (43 cm) | 18.5 in (46 cm) | | | |
| SEAT SURFACE HEIGHTS | 14 in (35 cm) | 15 in (51 cm) | 16 in (55 cm) | 17 in (59 cm) | 18 in (45 cm) | 19 in (48 cm) |
| WHEELS - MAG WHEELS | FRONT 20 in (35 cm) or 22 in (55 cm) | | | | | |
| | REAR 20 in (35 cm) or 22 in (55 cm) | | | | | |
| WHEELS - CASTER WHEELS | FRONT Aluminum Caster Wheels 5 in (12 cm) | | | | | |
| | FRONT Total Lock Caster Wheels 5 in (12 cm) | | | | | |
| | REAR Total Lock Caster Wheels 5 in (12 cm) | | | | | |
| LIFETIME WARRANTY FRAME | 16-Gauge Tubular Steel Frame 10 YEARS | | | | | |
| STANDARD PADDING | Upper & Lower Side Pads, Back Pad with Cover | | | | | |



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

> Click here to see accessories, color options, order forms, and more!

Available with solid seat pan option.



LIBRA CUSHION

Foam

HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.



MULTI CUSHION

Three-part foam combination Cooling gel foam and immersive support.



VI CUSHION

Ischial insert Supersoft viscoelastic foam insert for envelopment and protection of bony prominences.





PSP CUSHION

Polymer Insert Molded-in polymer gel provides extra skin protection and shear reduction.

brodaseating.com | 1-844-552-7632