

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy maneuverability and safe patient adjustments
- Reduces risk of injury and falls with supportive, stable positioning



Full tilt and recline functionality provides superior positioning



Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt-in-space enhances comfort



Comfort Tension Seating® increases support as well as sitting tolerance

Enhanced Positioning with refined features

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 38°						
DEGREES OF RECLINE	Up to 90° (Trendelenburg capable)						
WEIGHT CAPACITY	Up to 350 lbs (159 kg)						
SEAT WIDTHS	16 in (40 cm)	18 in (45 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)	28 in (71 cm)
SEAT DEPTHS	15.5 in (39 cm)	17 in (43 cm)	18.5 in (45 cm)	20 in (58 cm)			
SEAT SURFACE HEIGHTS (ADJUSTABLE)	17 in (43 cm)	18 in (45 cm)	19 in (48 cm)	20 in (58 cm)	21 in (53 cm)	22 in (55 cm)	
WHEELS - CASTER WHEELS	FRONT Total Lock Caster Wheels 5 in (12.5 cm)						
	REAR Fixed Multi-Directional Locking 8 in (20 cm)						
LIFETIME WARRANTY FRAME	16-Gauge Tubular Steel Frame 10 YEARS						
PRESSURE-REDUCING BACK SUPPORT	Integrated with Comfort Tension Seating®						



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)



Also available with rear-mounted mag wheel option