

Patient, Facility, & Caregiver Benefits

- *Increases quality of life and engagement with surroundings*
- *Allows for easy maneuverability and safe patient adjustments*
- *Reduces risk of injury and falls with supportive, stable positioning*




Full tilt and recline functionality provides superior positioning



Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt-in-space enhances comfort



Comfort Tension Seating® increases support as well as sitting tolerance



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)

Enhanced Positioning with refined features

SPECIFICATIONS

| | | | | | | | |
|-----------------------------------|---------------------------------------------------|---------------|-----------------|---------------|---------------|---------------|---------------|
| INFINITELY ADJUSTABLE TILT | Up to 38° | | | | | | |
| DEGREES OF RECLINE | Up to 90° (Trendelenburg capable) | | | | | | |
| WEIGHT CAPACITY | Up to 350 lbs (159 kg) | | | | | | |
| SEAT WIDTHS | 16 in (40 cm) | 18 in (45 cm) | 20 in (50 cm) | 22 in (55 cm) | 24 in (61 cm) | 26 in (66 cm) | 28 in (71 cm) |
| SEAT DEPTHS | 15.5 in (39 cm) | 17 in (43 cm) | 18.5 in (45 cm) | 20 in (50 cm) | 22 in (55 cm) | | |
| SEAT SURFACE HEIGHTS (ADJUSTABLE) | 17 in (43 cm) | 18 in (45 cm) | 19 in (48 cm) | 20 in (50 cm) | 21 in (53 cm) | 22 in (55 cm) | |
| WHEELS - CASTER WHEELS | FRONT Total Lock Caster Wheels 5 in (12.5 cm) | | | | | | |
| | REAR Fixed Multi-Directional Locking 8 in (20 cm) | | | | | | |
| LIFETIME WARRANTY FRAME | 16-Gauge Tubular Steel Frame 10 YEARS | | | | | | |
| PRESSURE-REDUCING BACK SUPPORT | Integrated with Comfort Tension Seating® | | | | | | |



Also available with rear-mounted mag wheel option