

Patient, Facility, & Caregiver Benefits

- Reduces the risk of patient and caregiver injury during bathing
- Removable arm supports make transfers safe and easy
- Adjustable components make for a comfortable bathing experience



Concave back option accommodates and supports gluteal shelf



Dynamic tilting and reclining permits safe and thorough bathing



Anterior and posterior tilt allows for easy user positioning



Rust-resistant frame and components are durable and long-lasting

Safety and Comfort for Hygiene Care

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 14° (with posterior tilt) Up to 8° (with anterior tilt)		
DEGREES OF RECLINE	Up to 65°		
WEIGHT CAPACITY	800 lbs. (362kg)		
SEAT WIDTHS	27 in (68 cm)	33 in (83 cm)	
SEAT DEPTHS	21 in (54 cm)	23 in (58 cm)	
SEAT SURFACE HEIGHT	24.5 in (63 cm)		
BACK HEIGHT	26 in (66 cm)		
LIFETIME WARRANTY FRAME	16-Gauge Tubular Steel Frame 10 YEARS		



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)



Need a lightweight option?

The Revive Shower Commode accommodates users up to 400 lbs.!