

Patient, Facility, & Caregiver Benefits

- Reduces the risk of patient and caregiver injury during bathing
- Removable arm supports make transfers safe and easy
- Adjustable components make for a comfortable bathing experience



Dynamic tilting and reclining permits safe and thorough bathing



The recline feature and removable arm supports make it easy for thorough bathing and skin assessments



Anterior and posterior tilt allows for easy user positioning



Rust-resistant frame and components are durable and long-lasting

Safety and Comfort for Hygiene Care

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 12° (with posterior tilt) Up to 10° (with anterior tilt)		
DEGREES OF RECLINE	Up to 67°		
WEIGHT CAPACITY	400 lbs. (181kg)		
SEAT WIDTHS	18 in (40 cm)	22 in (55 cm)	
SEAT DEPTHS	18 in (40 cm)		
SEAT SURFACE HEIGHT	24.5 in (63 cm)		
BACK HEIGHT	26 in (66 cm)		
LIFETIME WARRANTY FRAME	16-Gauge Tubular Steel Frame 10 YEARS		



STRAP COLOR OPTIONS



Click here to see
accessories, color options,
order forms, and more!



Need a bariatric option?

*The Revive Bariatric Shower Commode
accommodates users up to 800 lbs.!*