

### **Patient, Facility, & Caregiver Benefits**

- Reduces the risk of patient and caregiver injury during bathing
- Removable arm supports make transfers safe and easy
- Adjustable components make for a comfortable bathing experience



*Dynamic tilting and reclining permits safe and thorough bathing*



*The recline feature and removable arm supports make it easy for thorough bathing and skin assessments.*



*Anterior and posterior tilt allows for easy user positioning*



*Rust-resistant frame and components are durable and long-lasting*

### Safety and Comfort for Hygiene Care

#### SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 12° (with posterior tilt) Up to 10° (with anterior tilt)		
DEGREES OF RECLINE	Up to 67°		
WEIGHT CAPACITY	400 lbs. (181kg)		
SEAT WIDTHS	18 in (40 cm)	22 in (55 cm)	
SEAT DEPTHS	18 in (40 cm)		
SEAT SURFACE HEIGHT	24.5 in (63 cm)		
BACK HEIGHT	26 in (66 cm)		
LIFETIME WARRANTY FRAME	16-Gauge Tubular Steel Frame 10 YEARS		



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)



## Need a heavy-duty option?

*The Revive Bariatric Shower Commode accommodates users up to 800 lbs.!*