

Patient, Facility, & Caregiver Benefits

- *Infinitely adjustable tilt makes it quick and simple to reposition the user*
- *Designed for safety, stability, and reduced risk of injury*
- *Provides stimulation and contentment - excellent for memory care*



Front-pivot seat tilt maintains user comfort with smooth, effortless tilt



Adjustable arm-rests aid in user comfort and positioning



Adjustable seat surface height allows for proper foot propulsion



Comfort Tension Seating® improves postural support and sitting tolerance

Enhanced Mobility & Freedom

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 18°						
WEIGHT CAPACITY	350 lbs. (158 kg)						
SEAT WIDTHS	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)	28 in (71 cm)
SEAT DEPTHS	15.5 in (37.5cm)	17 in (43 cm)	18.5 in (47 cm)	20 in (51 cm)			
SEAT SURFACE HEIGHTS (ADJUSTABLE)	12.5 in (32 cm)	13.5 in (35 cm)	14.5 in (37 cm)	15.5 in (40 cm)	16.6 in (42 cm)	17.5 in (45 cm)	
WHEELS - MAG WHEELS	FRONT or REAR MOUNT 20 in (50 cm) or 22 in (55 cm)						
WHEELS - CASTER WHEELS	FRONT Aluminum Non-locking Caster Wheels						
	FRONT Fixed Caster Wheels 5 in (13cm)						
	REAR Total lock caster wheels 5 in (13cm)						
LIFETIME WARRANTY FRAME	16-Gauge tubular steel frame 10 YEARS						
PRESSURE-REDUCING SEAT & BACK SUPPORT	Integrated with Comfort Tension Seating®						



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)

