

## Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Promotes safe patient handling with enhanced positioning
- Proper seating alignment reduces patient fatigue and discomfort



*Tilt functionality provides additional pressure offloading*



*Front-pivot seat tilt maintains a forward line of sight, enhancing engagement.*



*Adjustable upper lateral supports fit a variety of users*



*Optional WC-19 Transportation Package enables safe patient travel*

*Shown with rear 14" mag wheel option*

### Positioning features with contemporary design

#### SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 40°					
INFINITELY ADJUSTABLE RECLINE	Up to 90° (Trendelenburg Capable)					
WEIGHT CAPACITY	350 lbs. (158 kg)					
SEAT WIDTHS	16 in (40 cm)	18 in (45 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)
SEAT DEPTHS	17 in (43cm)	18.5 in (46 cm)	20 in (50 cm)	22 in (56 cm)		
SEAT SURFACE HEIGHTS	19.5 in (48cm)	20.5 in (51 cm)	22 in (55 cm)	23.5 in (59 cm)		
WHEELS - MAG WHEELS	<b>REAR MOUNT</b> 14 in (35 cm), 20 in (50 cm)					
	<b>MID MOUNT</b> 14 in (35 cm)					
LIFETIME WARRANTY FRAME	16-Gauge tubular steel frame <b>10 YEARS</b>					
PRESSURE-REDUCING SEAT & BACK SUPPORT	Comfort Tension Seating® Back & Seat Supports Optional Matrix® Seating Series Cushions					
STANDARD PADDING	Neckrest, upper laterals, side pads, sole and calf pads, seat and back pads					



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)



## Synthesis Rehab

Check out the Synthesis Rehab which is Medicare reimbursable! use HCPCS E1161. (USA Only)