BRODA®

Synthesis Rehab Wheelchair

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Promotes safe patient handling with enhanced positioning
- Proper seating alignment reduces patient fatigue and discomfort



Shown with rear 20" mag wheels and matrix PS Cushion



Tilt and recline functionality provide pressure offloading.



Front-pivot seat tilt maintains a forward line of sight, enhancing engagement.



Adjustable upper lateral supports fit a variety of users.



Solid seat pan provides support for a variety of cushions.

BRODA



Dynamic positioning throughout rehabilitation

SPECIFICATIONS

INFINITELY ADJUSTABLE TILT	Up to 40°					
INFINITELY ADJUSTABLE RECLINE	Up to 90°					
WEIGHT CAPACITY	350 lbs. (158 kg)					
SEAT WIDTHS	16 in (40 cm)	18 in (40 cm)	20 in (50 cm)	22 in (55 cm)	24 in (61 cm)	26 in (66 cm)
SEAT DEPTHS	17 in (43cm)	18.5 in (46 cm)	20 in (50 cm)	22 in (56 cm)		
SEAT SURFACE HEIGHTS	19.5 in (48cm)	20.5 in (51 cm)	22 in (55 cm)	23.5 in (59 cm)		
WHEELS - MAG WHEELS	REAR MOUNT 14 in (35 cm), 20 in (50 cm)					
	MID MOUNT 14 in (35 cm)					
WHEELS - CASTER WHEELS	FRONT total lock or non-locking caster wheels 5 in (13cm)					
	REAR multi-directional locking caster wheels 6 in (15cm)					
LIFETIME WARRANTY FRAME	16-Gauge tubular steel frame 10 YEARS					
PRESSURE-REDUCING SEAT & BACK SUPPORT	Comfort Tension Seating [®] Back & Seat Supports Optional Matrx [®] Seating Series Cushions					



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

> Click here to see accessories, color options, order forms, and more!



LIBRA CUSHION

Foam

HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.



MULTI CUSHION

Three-part foam combination Cooling gel foam and immersive support.



VI CUSHION

Ischial insert Supersoft viscoelastic foam insert for envelopment and protection of bony prominences.





PSP CUSHION

Polymer Insert Molded-in polymer gel provides extra skin protection and shear reduction.

brodaseating.com | 1-844-552-7632