

## **Patient, Facility, & Caregiver Benefits**

- Provides necessary stimulation for conditions like dementia, Alzheimer's disease, Down syndrome, and Autism.
- Gliding motion promotes feelings of contentment and safety
- Promotes proper spinal posture and pelvic positioning




*Gentle gliding motion*



*Stimulates active patients and reduces the need to wander*



*Auto-locking feature increases safety and prevents falls*



*Comfort Tension Seating<sup>®</sup> increases support and sitting tolerance*

*Relaxing serenity with steady, calming motion*

## SPECIFICATIONS

<b>WEIGHT CAPACITY</b>	250 lbs* (114 kg)
<b>WHEELS</b>	Rear transport wheel, wheels for non-patient transport**
<b>LIFETIME WARRANTY FRAME</b>	16-Gauge Tubular Steel Frame 10 YEARS
<b>COMMON ACCESSORIES</b>	Padded armrests, full or half tray, adjustable armrests, contoured headrest

\*Adjustments available to accommodate larger patients

\*\*The Tranquille Glider is not a transport chair

MODEL	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	BACK HEIGHT
SMALL LOW-BACK 100-10AL	18 in (45 cm)	17 in (43 cm)	16 in (40 cm)	26 in (66 cm)
LOW-BACK 100-15AL	20 in (50 cm)	17 in (43 cm)	18 in (45 cm)	26 in (66 cm)
HIGH-BACK 100-20AL	20 in (50 cm)	17 in (43 cm)	18 in (45 cm)	32 in (81 cm)



Additional accessories are available, including arm supports, trays, positioning straps, and more!

[Click here to see accessories, color options, order forms, and more!](#)

