

## **Patient, Facility, & Caregiver Benefits**

- Pressure redistribution that aids in maintaining skin integrity
- Improved air flow, decreasing heat and moisture build up
- Accommodation for conditions such as kyphosis and lordosis



*Shown with the Acute Care Configuration (ACC) full vinyl wrap option*



*Anterior and posterior tilt aids in safe patient handling*



*Swing Away, Removable Arm Supports facilitate side-transfers*



*Manual foot pump allows safe, no lift tilt positioning*



*Concave back for gluteal shelf improves user support & dignity*

### Bariatric wheelchair care redefined

#### SPECIFICATIONS

<b>INFINITELY ADJUSTABLE TILT</b>	Up to 10° Anterior, Up to 30° Posterior					
<b>DEGREES OF RECLINE</b>	Up to 45°					
<b>WEIGHT CAPACITY</b>	Up to 600 lbs (272 kg)					
<b>SEAT WIDTHS</b>	24 in (61 cm)	26 in (66 cm)	28 in (71 cm)	30 in (76 cm)		
<b>SEAT DEPTHS</b>	21 in (53 cm)	23 in (58 cm)				
<b>SEAT SURFACE HEIGHTS</b>	20 in (50 cm)	21 in (53 cm)	22 in (55 cm)	23 in (58 cm)	24 in (60 cm)	
<b>WHEELS - CASTER WHEELS</b>	<b>FRONT</b> Total Lock Caster Wheels 6 in (15 cm)					
	<b>REAR</b> Fixed Multi-Directional Locking 8 in (20 cm)					
<b>LIFETIME WARRANTY FRAME</b>	16-Gauge Tubular Steel Frame 10 YEARS					
<b>PRESSURE-REDUCING BACK SUPPORT</b>	Integrated with Comfort Tension Seating*					
<b>BACK SUPPORT</b>	Upper & Lower Side Pads, Back Pad with Cover					



Additional accessories are available, including foot supports, trays, positioning straps, padding packages, and more!

[Click here to see accessories, color options, order forms, and more!](#)

